

# Ligators

MACQUARIE MEDICAL SYSTEMS  
CATALOGUE VERSION 11



Ligator Bands	3
Ligator, Haemorrhoidal	3
McGiveny Ligator Forcep	3



# FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.


	<p><b>THEATRE (GERMANY)</b></p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
	<p><b>GERMANY</b></p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
	<p><b>SUPERIOR</b></p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
	<p><b>WARD</b></p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



### MCGIVENY LIGATOR FORCEP

McGiveny Ligator Forceps are used to grasp the ligator bands used in haemorrhoid removal.

	HIPP	KLINI	ARMO
19cm			A3605



### LIGATOR, HAEMORRHOIDAL

Rudd Haemorrhoid Ligator is an instrument used in one process for the removal of haemorrhoids. It uses a thumb plunger system to release the bands.

	HIPP	KLINI	ARMO
With cone		K33.350.31	



### LIGATOR, HAEMORRHOIDAL

Rudd Haemorrhoid Ligator is an instrument used in one process for the removal of haemorrhoids. It uses a thumb plunger system to release the bands.

	HIPP	KLINI	ARMO
With cone			A1990



### LIGATOR BANDS

	HIPP	KLINI	ARMO
NON-STERILE			
Vial of 100			A1991



301 Catherine St,  
Leichhardt NSW 2040

Phone: (02) 9692 7921  
[www.macquariemed.com.au](http://www.macquariemed.com.au)