

Biopsy Instruments

MACQUARIE MEDICAL SYSTEMS
CATALOGUE VERSION 11



Biopsy Forceps

Tischler 3
Yoemann 3

Biopsy Punches

Keyes 3

Biopsy Punches, Single Use

Accuderm Punches 4



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

 <p>ANTON HIPPI Instruments & Implants</p>	<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
 <p>KLINI™ instruments by anton hipp GERMANY</p>	<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
 <p>ARMO surgical instruments</p>	<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
 <p>Sayco first aid instruments</p>	<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>

Biopsy Forceps



TISCHLER

Tischler Biopsy Forceps are obstetric/gynaecological instruments most commonly used to obtain tissue samples.

	HIPP	KLINI	ARMO
STANDARD (7 X 3mm JAW)			
25cm shaft	2.550.07		A6030
BABY (4.2 X 2.3mm JAW)			
21cm shaft	2.550.04		A6031



YOEMANN

Yoemann Biopsy Forceps are pistol grip forceps used for obtaining rectal biopsies.

	HIPP	KLINI	ARMO
35cm			A6032
42cm			A6033

Biopsy Punches



○
2mm
Diametre

KEYES

Keyes Biopsy Punch is a dermal diagnostic instrument used to remove skin samples for biopsy.

	HIPP	KLINI	ARMO
2mm Ø			
10.5cm			A6034



Biopsy Punches, Single Use

STERILE ACCUDERM PUNCHES

Acu-Punch® is made in the U.S.A. — and the ONLY punch available in a complete range of sizes, from 2mm to 12mm. Acu-Punch® is always razor sharp and ready to use. The handle is ribbed for maximum control, even when wearing gloves. Each punch comes individually wrapped in an easy peel pouch. 50 punches per box.

Acu-Punch®	
50 PUNCHES PER BOX	
2mm	P250
3mm	P350
4mm	P450
5mm	P550
6mm	P650
8mm	P850
10mm	P1050
12mm	P1250

Instrument quality level



Anton Hipp
Theatre (Germany)

Klini
Germany

Armo
Superior

Sayco
Ward





301 Catherine St,
Leichhardt NSW 2040

Phone: (02) 9692 7921
www.macquariemed.com.au