

Hollowware Polypropylene

MACQUARIE MEDICAL SYSTEMS
CATALOGUE VERSION 11



Medicine Measures

Graduated Measure	3
Lid	3

Kidney Dishes

Heavy Duty	3
Standard	3

Instrument Trays

Heavy Duty	4
------------	---

Bowls

Heavy Duty	4
Standard	4



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

	<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
	<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
	<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
	<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>

Medicine Measures

GRADUATED MEASURE

These washable graduated measures are clear in colour and come in 2 sizes with capacities of 30ml and 60ml.

	Code
CLEAR, PACK OF 100	
2.5ml - 30ml	PMM30
CLEAR, PACK OF 25	
5ml - 60ml	PMM60

LID

The Lids are available to fit both the 30ml and 60ml measure and sold in a pack of 100.

	Code
CLEAR, PACK OF 100	
To suit 2.5ml-30ml & 5ml-60ml	PMM30L

Kidney Dishes

HEAVY DUTY

Available in 700 & 1000ml.

	Code
GREEN	
700ml	DK700GU
1000ml	DK1000GU

STANDARD

Available in 700ml dishes and packs of 25.

	Code
GREEN	
700ml	DK700GS





Instrument Trays

HEAVY DUTY

These instrument trays are heavy duty for repeated use. They are made in Australia and are available in both solid and perforated.

	CODE
GREEN SOLID	
220 x 140 x 20mm	CBF2214-2
245 x 195 x 40mm	CBF2419-4
GREEN PERFORATED	
180 x 150 x 30mm	CBF1815-3
270 x 150 x 30mm	CBF2715-3

Bowls



HEAVY DUTY

These autoclavable bowls are designed for repeated use, have permanent graduated markings and are Australian made. They range in size from 150ml to 1500ml.

	CODE
GREEN	
150ml (100 x 50mm)	PB100GU
500ml (140 x 70mm)	PB140GU
1000ml (185 x 90mm)	PB185GU
2500ml (240 x 110mm)	PB240GU



STANDARD

These washable bowls are for single patient use, they have permanent graduated measurement markings and are Australian made. They range in size from 250ml to 2000ml. There is also a denture cup available in this range and they are all available in packs of 25.

	CODE
GREEN, PACK OF 25	
250ml (95 x 50mm)	PB250G
500ml (115 x 60mm)	PB500G
1000ml (145 x 80mm)	PB1000G
2000ml (195 x 90mm)	PB2000G

Instrument quality level



Theatre (Germany)

Germany

Superior

Ward



Hollowware Polypropylene

05



301 Catherine St,
Leichhardt NSW 2040

Phone: (02) 9692 7921
www.macquariemed.com.au