

## Sphygmomanometers

MACQUARIE MEDICAL SYSTEMS  
CATALOGUE VERSION 11



### Aneroid

Econom - Economy Aneroid	2
Home Blood Pressure Kit	2

### Digital

Rossmax BP Automatic Upper Arm Bluetooth Monitor	3
Rossmax BP Automatic Upper Arm Monitor (Fuzzy Logic)	3
Rossmax BP Automatic Wrist Monitor	2
Rossmax BP Large Cuff 34-46cm	3

### Palm Held

Palm- HandHeld Aneroid	4
Palm Lite - Hand Held Aneroid	4

### Regal

Desk - Aneroid Sphygmomanometer Regal	4
Mobile - Aneroid Sphygmomanometer Regal	5
Wall Mounted - Aneroid Sphygmomanometer Regal	5



## Aneroid

### **ECONOM - ECONOMY ANEROID**

- Economically priced precision sphygmomanometer
- 300mmHg manometer with no pin stop mechanism.
  - Includes Adult size nylon fabric cuff with velcro.
  - Complete with heavy duty 2 tube inflation bladder, bulb and valve.
  - Supplied with zippered leatherette storage case.

	<b>ABN</b>
Kit	<b>SGAS-051</b>

### **HOME BLOOD PRESSURE KIT**

Precision sphygmomanometer for professional and home use

- Designed for professional clinical use or for easy home self monitoring by individuals.
- Comes with Nylon D-ring cuff and attached stethoscope.
- D-ring cuff fits arms 25 to 35cm in circumference.
- 300mmHg manometer with no pin stop mechanism.
- Complete with heavy duty 2 tube inflation bladder, bulb and valve.
- Supplied with zippered leatherette storage case.

	<b>ABN</b>
Kit	<b>SGAS-061</b>

## Digital

### **ROSSMAX BP AUTOMATIC WRIST MONITOR**

Features:

- Smart Sense technology
- Irregular Heartbeat Detection (IHB)
- Hypertension risk indication
- Average of last 3 readings
- 120 Memories for 2-users with date and time indication
- Movement detection
- Jumbo size screen
- PC link
- Lifetime calibration
- One-touch automatic operation

	<b>ROSSMAX</b>
Kit	<b>RMBQ75</b>

## Digital



### ROSSMAX BP AUTOMATIC UPPER ARM MONITOR (FUZZY LOGIC)

Automatic BP monitor with lifetime calibration.

- Real Fuzzy Logic Technology
- Cuff Detection
- Movement Detection
- Hypertension Risk Indication
- Irregular Heartbeat Detection
- Average of Last 3 readings
- 120 Memories for 2 users with date and time
- Large LCD
- One Touch Operation

#### ROSSMAX

Kit	<b>RMX3</b>
-----	-------------



### ROSSMAX BP AUTOMATIC UPPER ARM BLUETOOTH MONITOR

Automatic BP Monitor with Arrhythmia detection and lifetime calibration.

- Bluetooth low energy technology
- Arrhythmia Detection (ARR)
- Atrial Fibrillation Detection
- Premature Contraction Detection
- Movement Detection
- Cuff Detection
- Hypertension Risk Indication
- Average of last 3 readings
- 7 day AM/PM Averaging
- Patented Latex Free Universal Cone Cuff 24-40cm

#### ROSSMAX

Kit	<b>RMX5BT</b>
-----	---------------



### ROSSMAX BP LARGE CUFF 34-46CM

No more annoying rotating and sliding when wrapping the cuff! The ergonomic, cone-shaped design of Rossmax Cone cuffs promotes blood pressure measurement to a whole new easy and delightful experience! This cuff is Latex free and Easy to apply and is available in small, medium and large.

#### ROSSMAX

Kit	<b>BPCF001</b>
-----	----------------



## Palm Held

### **PALM- HANDHELD ANEROID**

- Deluxe handheld model sphygmomanometer
- Deluxe Sphygmomanometer palm style with 300mmHg with no pin stop mechanism.
  - Includes Adult size nylon fabric cuff with velcro.
  - Complete with heavy duty 1 tube inflation bladder.
  - Chrome plated inflation valve.
  - Unit supplied with deluxe zippered nylon storage case.

	<b>ABN</b>
Kit	<b>SGAS-071LT</b>



### **PALM LITE - HAND HELD ANEROID**

- Handheld sphygmomanometer with precision performance and a modern lightweight design
- Precision crafted 300mmHg palm style manometer with enhanced shock resistance housing.
  - Quick release air valve button.
  - Quick release luer slip connector.
  - Includes Adult size nylon fabric cuff with velcro.
  - Complete with heavy duty 1 tube inflation bladder.
  - Unit supplied with deluxe zippered nylon storage case.

	<b>ABN</b>
Kit	<b>SGAS-071</b>



## Regal

### **DESK - ANEROID SPHYGMOMANOMETER REGAL**

- Portable clock type Sphygmomanometer for professionals.
- Large 150mm faceplate 300mmHg manometer with no pin stop mechanism.
  - Easy to read black background with white writing.
  - Special non-reflective manometer glass.
  - Extra strong reinforced base design.
  - Includes Adult size nylon fabric cuff with velcro.
  - Complete with heavy duty 2 tube inflation bladder, bulb and valve.
  - Convenient under manometer cuff storage.
  - 120cm extendable coiled tubing included.

	<b>REGAL</b>
Kit	<b>SGAS-081D</b>

## Regal



### MOBILE - ANEROID SPHYGMOMANOMETER REGAL

- Clock type Sphygmomanometer on mobile Floor Stand.
- Large 150mm faceplate 300mmHg manometer with no pin stop mechanism.
  - Easy to read black background with white writing.
  - Special non-reflective manometer glass.
  - Two-piece height adjustable stand.
  - Five caster wheel base ensures easy movement and high stability.
  - Built-in cuff storage.
  - Stand is finished in durable chrome and clinical white enamel.
  - Includes Adult size nylon fabric cuff with velcro.
  - Complete with heavy duty 2 tube inflation bladder, bulb and valve.
  - 240cm extendable coiled tubing included.
  - All necessary assembly tools included.

---

**REGAL**

Kit **SGAS-081M**

---



### WALL MOUNTED - ANEROID SPHYGMOMANOMETER REGAL

- Portable clock type Sphygmomanometer for professionals.
- Alternative presentation of the CLOCK aneroid sphygmomanometer.
  - Attaches to 180° swivel wall mounting bracket.
  - Large 150mm faceplate 300mmHg manometer with no pin stop mechanism.
  - Easy to read black background with white writing.
  - Special non-reflective manometer glass.
  - Extra strong reinforced storage basket.
  - Includes Adult size nylon fabric cuff with velcro.
  - Complete with heavy duty 2 tube inflation bladder, bulb and valve.
  - 240cm extendable coiled tubing included.

---

**REGAL**

Kit **SGAS-081W**

---



301 Catherine St,  
Leichhardt NSW 2040

Phone: (02) 9692 7921  
[www.macquariemed.com.au](http://www.macquariemed.com.au)