

## Thermometers

MACQUARIE MEDICAL SYSTEMS  
CATALOGUE VERSION 11



Flexible Tip Thermometer & Probe	3
Infrared Ear Thermometer	2
Non-Contact Temple Thermometer	3
Probe Cover	2



°C Reading

°F Reading



# Infrared Ear Thermometer

## INFRARED EAR THERMOMETER

ROSSMAX Infrared Ear Thermometer features

- 1 second measurement or continuous scanning modes
- Fever alarm
- 9 memories
- °C / °F switchable display readout with backlight
- Auto shut-off
- Low battery indicator
- High accuracy to ASTM Clinical Standard
- Self diagnosis when device is abnormal

	<b>ROSSMAX</b>
Kit	<b>RMRA600</b>

## PROBE COVER

ROSSMAX Infrared Ear Thermometer probe cover features:

- Patented PLUG & SECURE probe covers
- Provide sanitary environment for patient temperature reading
- Sold in packs of 25

	<b>ROSSMAX</b>
Packs of 25	<b>RMRA601</b>

## Non-Contact Temple Thermometer



### NON-CONTACT TEMPLE THERMOMETER

ROSSMAX Non Contact Thermometer.

- 1 second measurement
- Fever alarm
- 9 memories
- °C / °F switchable display readout with backlight
- Auto shut-off
- Low battery indicator
- High accuracy to ASTM Clinical Standard
- Self diagnosis when device is abnormal

	<b>ROSSMAX</b>
Kit	<b>RMHA500</b>

## Flexible Tip Thermometer & Probe



### FLEXIBLE TIP THERMOMETER & PROBE

Thermometer Features

- Accurate measurement in 10 seconds
- Fever alarm
- Waterproof
- Auto shut off
- °C / °F switchable display readout
- High accuracy to ASTM Clinical Standard
- Low battery indicator
- Self diagnosis when device is abnormal

	<b>ROSSMAX</b>
Kit	<b>RMTG100</b>



301 Catherine St,  
Leichhardt NSW 2040

Phone: (02) 9692 7921  
[www.macquariemed.com.au](http://www.macquariemed.com.au)