


 **Macquarie**
Medical Systems

gmmedical
Product Catalogue




GM Medical's extensive range of equipment continues the tradition of industry leadership.

With the recent inclusion into the ParagonCare Group GM Medical has expanded the range of equipment to offer significantly more products into the healthcare environment.

GM Medical has developed a comprehensive equipment range that combines advanced functionality and contemporary design with quality and reliability.


Using a variety of finishes and colour combinations the product range offers unlimited applications throughout the healthcare environment, from the stringent infection control requirements of the operating theatre and CSSD right through to the Aged Care facilities and Consulting Rooms with colour coordinated environments.

Instrument Trolleys	01
Equipment Trolleys	04
Spectra Trolleys	08
Sixty Forty Materials Management	09
CSSD Equipment	10
Operating Theatre	12
General Ward	16
Maternity	24
Linen and Supply	26
Food Services	31
Consulting Rooms	32
Lighting	37
Product Guide	38



GM Medical instrument trolley design methods are class leaders in providing high quality carts for the surgical and clinical environments.

With a focus on durability features such as double folded safety edges, under shelf bracing and castor mount inserts as standard. These ensure safety and stability at all times.



Shelf Rails
12mm stainless steel rails
Fully welded stainless steel shelf
Double safety folds to all edges

Under Shelf Bracing
Full under bracing providing maximum stability and strength

Drawer Runners
Full extension, easy glide metal runners with drawer retention system

Stainless Steel Construction
Satin finish - 304 grade stainless steel throughout

Shelf - Stainless Steel
Double safety folds to all edges

Castors
75mm ball bearing fork castor
Soft grey rubber wheel
Expanding adaptor for easy removal

NO RAIL



AX 001
Instrument Trolley
490 x 490 x 900mm



AX 003
Instrument Trolley
600 x 490 x 900mm



AX 006
Instrument Trolley
750 x 490 x 900mm



AX 009
Instrument Trolley
900 x 490 x 900mm
Other models available
AX 0015
1000 x 600 x 900mm
AX 0018
1200 x 600 x 900mm

RAIL 3 SIDES



AX 033
Instrument Trolley
490 x 490 x 900mm



AX 036
Instrument Trolley
600 x 490 x 900mm



AX 039
Instrument Trolley
750 x 490 x 900mm



AX 042
Instrument Trolley
900 x 490 x 900mm
Other models available
AX 0048
1000 x 600 x 900mm
AX 0051
1200 x 600 x 900mm

ONE DRAWER



AX 054
Instrument Trolley
490 x 490 x 900mm



AX 057
Instrument Trolley
600 x 490 x 900mm



AX 060
Instrument Trolley
750 x 490 x 900mm



AX 063
Instrument Trolley
900 x 490 x 900mm

TWO DRAWER



AX 066
Instrument Trolley
490 x 490 x 900mm



AX 069
Instrument Trolley
600 x 490 x 900mm



AX 072
Instrument Trolley
750 x 490 x 900mm



AX 075
Instrument Trolley
900 x 490 x 900mm



AX 076
Instrument Trolley
490 x 490 x 900mm
Stainless steel construction
3 Drawers - 125mm deep
125mm castors

Other models available
AX 077
600 x 490 x 900mm
3 Drawers - 125mm deep
AX 0780
750 x 490 x 900mm
3 Drawers - 125mm deep
AX 079
900 x 490 x 900mm
3 Drawers - 125mm deep



AX 080
Instrument Trolley
490 x 490 x 900mm
Stainless steel construction
4 Drawers - 125mm deep
125mm castors


Other models available
AX 081
600 x 490 x 900mm
4 Drawers - 125mm deep
AX 082
750 x 490 x 900mm
4 Drawers - 125mm deep
AX 083
900 x 490 x 900mm
4 Drawers - 125mm deep



AX 073
Instrument Trolley
900 x 490 x 900mm
Stainless steel construction
2 half drawers - 125mm deep
75mm castors



AX 107
Resuscitation Trolley
750 x 490 x 900mm
Stainless construction
4 half drawers - 125mm deep
125mm castors



The range of Equipment Trolleys is extensive and all trolleys are designed featuring a solid tubular frame which ensures rigidity at all times. This allows the drawers to be mounted securely and to operate efficiently. All drawer runner systems are full ball bearing extension type featuring a push lock mechanism to prevent the drawer from accidentally opening during transport. As with the instrument trolleys the

shelves feature double fold and under shelf bracing for a stable work platform. Our range of Equipment trolleys are available in all stainless steel or can be colour coordinated with powder coated drawers and panels. All trolleys are fitted with 125mm precision bearing heavy duty castors with optional corner buffers available.



AX 108
Resuscitation Trolley - Standard
750 x 490 x 900mm
Stainless construction
4 Half drawers - 125mm deep
1 Full width drawer - 125mm deep
125mm castors



AX 110
Resuscitation Trolley
750 x 490 x 900mm
Stainless construction
6 Half drawers - 125mm deep
1 Full width drawer - 250mm deep
125mm castors



AX 078
Multi Drawer Trolley
900 x 490 x 900
Stainless steel construction
4 Drawers - 100mm deep
2 Drawers - 150mm deep
125mm castors



AX 104
Anaesthetic Trolley
600 x 490 x 1000mm
Stainless construction
2 Drawers - 125mm deep
2 Drawers - 250mm deep
125mm castors



AX 102
Anaesthetic Trolley
600 x 490 x 900mm
Stainless construction
3 Drawers - 125mm deep
1 Drawer - 250mm deep
125mm castors



AX 105
Anaesthetic Trolley
900 x 490 x 900mm
Stainless construction
6 Full width drawers
125mm castors

AX 114 (not pictured)
Anaesthetic Trolley
600 x 490 x 900mm
Stainless construction
5 Full width drawers -
125mm deep
125mm castors



AX 117
Emergency Equipment Trolley
 750 x 490 x 900mm
 Stainless steel frame
 Powder coated panels and drawers
 4 Drawers - 105mm deep
 1 Drawer - 225mm deep
 125mm castors



AX 115
Emergency Equipment Trolley
 750 x 490 x 900mm
 Stainless steel frame
 Powder coated panels and drawers
 5 Drawers - 125mm deep
 125mm castors



AX 145
IV Set-up Trolley
 900 x 490 x 1122mm
 5 Drawers - 125mm deep
 2 Half drawers - 250mm deep
 125mm castors



AX 119
Emergency Equipment Trolley
 750 x 490 x 1120mm
 Stainless steel frame
 Powder coated panels and drawers
 3 Drawers - 125mm deep
 2 Drawers - 250mm deep
 125mm castors

MEDICATION TROLLEYS



AX 676
Medication Trolley - 3 Drawer
 Stainless construction
 600 x 490 x 1100mm
 Lift up shelf
 3 Drawers - 275mm deep
 125mm castors



AX 677
Medication Trolley - 1 Drawer
 Powder coated finish
 600 x 490 x 900mm
 1 Webster drawer
 275mm deep - lockable
 Lift up shelf
 75mm castors

TROLLEY ACCESSORIES




AX 120
Drawer insert containers - 12 units



AX 738
Self closing drawer runner



AX 132
Storage panel - 4 sml / 3 med / 1 lrg



AX 126
Drop down work surface



AX 123
Vertical oxy bottle holder and clamp



AX 129
IV pole and clamp



AX 135
4 x pin power board




AX 143
Push handle
 to suit 400mm wide trolley



AX 142
Waste bin and bracket

SPECTRA TROLLEYS



The Spectra range is a trolley system that combines mobility and drawer flexibility into a lightweight colour coded package that has applications throughout all healthcare environments.

The unique drawer systems is based on a removable 60 x 40 tray system that offers 3 different depths and the option of integrated drawer dividers.

This allows flexibility with the trolley applications, particularly in the ability to restock the carts with transfer tray systems.

The shelves are one piece moulded polymer which are easy to clean and eliminates infection control issues due to the rounded internal mould design, whilst integrating a raised lip to three sides to prevent items falling off.



AX S0001
2 Shelf Trolley
700 x 500 x 1000mm
Work surface area 620 x 460mm
125mm castors

Push handle optional on all



AX S0010
1 Drawer Trolley
700 x 500 x 1000mm
AX M0601 - 1 x 100mm deep tray
Interlock dividers optional
Work surface area 620 x 460mm
125mm castors



AX S0020
2 Drawer Trolley
700 x 500 x 1000mm
AX M0601 - 2 x 100mm deep trays
Interlock dividers optional
Work surface area 620 x 460mm
125mm castors



AX S0030
Multi Drawer Trolley
700 x 500 x 1000mm
AX M0605 - 1 x 50mm deep tray
AX M0601 - 2 x 100mm deep trays
AX M0602 - 2 x 200mm deep trays
Interlock dividers optional
Work surface area 620 x 460mm
125mm castors

COLOUR OPTIONS



60 x 40 Tray System



50mm Trays
AX M0605 50mm deep tray
AX M1665 Interlock divider 600mm
AX M1645 Interlock divider 400mm



100mm Trays
AX M0601 100mm deep tray
AX M1661 Interlock divider 600mm
AX M1641 Interlock divider 400mm



200mm Trays
AX M0602 200mm deep tray
AX M1662 Interlock divider 600mm
AX M1642 Interlock divider 400mm

SIXTY | FORTY



The 60 x 40 materials management system is recognised throughout the world as the industry leader in hospital and healthcare storage. GM Medical brings these advanced storage solutions to Australia.

The 60 x 40 concept is based on a tray size of 600 x 400mm and a choice of tray depths of 50, 100, 200mm. This system allows versatility and flexible

configuration, as the trays are interchangeable with the range of GM cabinets, shelving systems and carts.

Choose your Configuration All GM Medical 60 x 40 trolleys provide you the opportunity to configure them to your requirements. Mix and match using our suggestions pictured as your starting point.



AX M2000
Mobile Storage Unit
900 x 650 1000mm
2 x 50 / 6 x 100 / 2 x 200mm deep trays
125mm castors - 2 brake
Push handle



AX M1400
Mobile Storage Unit
665 x 460 x 975mm
4 x 100 / 1 x 200mm deep trays
100mm castors - 2 brake
Push handle



AX M2500
Mobile Storage Unit
Stainless Steel Construction
600 x 460 x 900mm
2 x 120mm deep drawers
1 x 200mm deep tray
125mm castors



AX M8000
Mobile Storage Unit
680 x 450 x 1950mm
1 x 50 / 4 x 100mm deep trays
2 x 200mm deep trays
125mm castors - 2 brake



AX M8500
Static Storage Unit
680 x 450 x 1950mm
1 x 50 / 4 x 100mm deep trays
2 x 200mm deep trays
Adjustable feet



Tray Stop System
Prevents accidental removal of trays.
Multi position holes enable flexible configuration.

CSSD EQUIPMENT

Our range of CSSD workstations and associated equipment are designed with flexibility in mind.

A large range of accessories, options and sizes are available and allows the workstation to be configured to meet a multitude of applications, optimising workflow, room layout and departmental efficiency.

The simple push button electric height adjustment ensures that the optimum working height is achieved to ensure correct operator's posture, providing comfort, and minimising occupational health and safety issues.

The workstations are available in six standard sizes, however should there be a requirement for

a particular size, we are able to custom manufacture.

Construction of all stainless steel, featuring a reinforced heavy duty top surface, ensures infection control issues are minimised. Stable base design ensures stability and provides clear access for operators legs and feet.

PICTURED WITH ACCESSORIES



Accessory Shelf
Constructed of chrome plated wire. 175mm wide. May be positioned at different heights if required.
Option - Shelf depth 300mm
Additional shelves may be added
AX 0672 1800 x 175mm
AX 0677 1800 x 300mm
AX 0673 1500 x 175mm
AX 0678 1500 x 300mm
AX 0674 1200 x 175mm
AX 0679 1200 x 300mm

AX 0685
Sterile Wrap Rack
This rack system can be fitted either side of the workstation and is removable and telescopic.
Option - The wrap holders have the option of being telescopic to accommodate wider rolls of sterile wrap.

AX 0687
Keyboard Shelf
Pull out keyboard shelf on quality roller drawer runners. May be positioned anywhere along the length of the workstation.

AX 0686
CPU Holder
Robust construction. Can be positioned anywhere beneath the workstation surface. Designed to for standard CPU units.
Option - custom size CPU holder

Additional Accessories (not pictured)
AX 0689 Lower Storage Shelf
AX 0690 Power Board Outlet 6 plug
AX 0691 Steri-Peel Holder

Electric Height Controls
Simple push button controls conveniently located at the corner of the workstation. Dual column electric height adjustment. Standard height range 800-1100mm. The height range is 300mm.
Option - Battery back up system

Size Options
AX 0666 1800 x 760mm
AX 0669 1800 x 900mm
AX 0667 1500 x 760mm
AX 0670 1500 x 900mm
AX 0678 1200 x 760mm
AX 0671 1200 x 900mm

BASIC UNIT



CSSD EQUIPMENT



AX 692
Case Cart
900 x 600 x 1040mm
Stainless steel construction
Fold out hinged doors to 270°
Internal adjustable wire shelves
Top rail push handle incorporated
Full bumper surround
200mm castors - 2 brake



AX 278
List Cart
1200 x 600 x 1430mm
Stainless steel construction
2 hinged doors / internal adjustable shelves
200mm castors / full bumper surround
In addition to the standard sizes offered, custom sizes are available for both List and Case Carts.

Cart Features



Door System
The doors swing open to a full 270° and ensures full unobstructed access. Flush fitting, heavy duty door latches are fitted. The doors have bracing rods welded to the door perimeter to ensure rigidity and durability. The door hinges are all stainless steel, powder coated white.
Shelf System
Each unit is supplied with an adjustable shelf system that allows the stainless steel wire shelf to be positioned as required. Additional shelves are available.



AX 0688
Preparation Table Adjustable Height
1800 x 900mm / height range 830 - 1130mm
Electrically operated height adjustment
Stainless steel work surface
Storage shelves / battery back up
100mm braking castors




AX 0684
Transportation Trolley
1200 x 600 x 900mm
Stainless steel construction
3 shelves with rails to three sides
125mm castor



AX 0681
Kinguard Drying Rack
1200 x 500 x 1500mm
Stainless steel construction
6 rails / 125mm castors



AX 0680
Packing Table
Stainless steel construction
1 drawer - 125mm deep
Pull out keyboard shelf



The GM Medical range of operating theatre equipment has a strong emphasis on being manufactured from stainless steel minimising infection control issues. All products in this range are fully welded and highly polished, therefore easy to clean and maintain.

This diverse range of products covers all clinical and surgical equipment requirements.



AX 106
Airway Trolley
750 x 490 x 900mm
Stainless construction
1 Drawer - 125mm deep
Clear Acrylic top



AX 4990
Plaster Trolley
900 x 490 x 900mm
Stainless construction
2 Drawers - 125mm deep
Drop down shelf / bucket hook



AX 252
Mayo Table
Height 820 - 1120mm
Tray size 560 x 410mm
Stainless construction
2 Leg flat base
75mm castors



AX 253
Mayo Table
Height 820 - 1120mm
Tray size 560 x 410mm
Stainless construction
2 Leg slide base
75mm castors



AX 254
Mayo Table
Height 820 - 1120mm
Tray size 560 x 410mm
Stainless construction
3 Leg base
75mm castors



AX 256
Mayo Table
Height 820 - 1120mm
Tray size 560 x 410mm
Stainless construction
4 Leg base
50mm castors



AX 021
Nesting Tables - Set of Three
AX 024
Nesting Table - Large
1000 x 600 x 900mm
AX 027
Nesting Table - Medium
940 x 510 x 880mm
AX 030
Nesting Table - Small
840 x 490 x 840mm



AX 111
Arm Table - Fixed Height
750 x 490 x 900mm
Stainless steel construction
Safety folds to all edges
Reinforced shelf
2 x 75mm castors



AX 113
Arm Table - Adjustable Height
750 x 490mm
Height range 740 - 1140
Stainless steel construction
Reinforced shelf
2 x 75mm castors
Radiolucent top



AX 260
Lithotomy Table
750 x 490mm Adjustable height
4 - 75mm castors



AX 273
Surgeons Stool



AX 277
Theatre Stool

AX 273

Surgeons Stool - Hydraulic

- Special operators stool for surgery
- Height adjustment by maintenance free hydraulic pump lift
- Safety base made of stainless steel with 5 antistatic twin castors
- Chromium plate foot ring
- Stroke 190mm
- Seat height range 530 - 720
- Back rest adjustment range 70mm
- Upper part swivelling & lockable
- Comfortable antistatic black upholstery
- Pair of upholstered OP - arm supports
- Base locking device

AX 277

Theatre Stool

- Height adjustment by safety gas spring controlled by stainless steel wobble plate
- Multi purpose stool for surgery, anaesthesia and outpatient department
- Safety base made of stainless steel with 5 twin castors for hard floors
- Round seat with antistatic black upholstery
- Seat height range 500 - 680mm
- Base locking device all round operable by step disc above base



AX 262
Surgeon Stool Hand Operated
Height 530 - 710mm
No backrest



AX 265
Surgeon Stool Hand Operated
Height 530 - 710mm
Backrest



AX 267
Surgeon Stool Foot Operated
Height 530 - 710mm
No backrest



AX 270
Surgeon Stool Foot Operated
Height range 530 - 710
Backrest



AX 2710
Saddle Stool
Hand operated height adjustment
Height range 530 - 730mm



AX 391
One Step Stool with hand rail
Aluminium tread plate
700 x 300 x 200mm
Handle height; 1000mm



AX 387
One Step Stool Aluminium Tread
400 x 300 x 200mm
Stainless steel frame
Aluminium tread plate
AX 390
One Step Stool Aluminium Tread
700 x 300 x 200



AX 405
Two Step Stool Aluminium Tread
Step Height: 200 / 400mm
Stainless steel frame
Aluminium tread plate



AX 214
Double Bowl Stand
Height 820mm
300mm (dia) Stainless steel bowl



AX 211
Single Bowl Stand
Height 820mm
300mm (dia) Stainless steel bowl



AX 248
Single Bowl Stand
Height 515mm
300mm (dia) Stainless steel bowl



AX 240
Kick Bucket Stand
Bucket and stand combination
Stainless steel construction
Separate components
AX 241 Kick Bucket Stand
AX 244 Kick Bucket - 11.4 litre
AX 247 Kick Bucket Lid



AX 701



AX 702

AX 701

ICU Work Station

Writing surface: 750 x 600mm
Height range: 700 - 1000mm
Powder coated
Tilting writing surface

AX 702

ICU Work Station

1500 x 600 x 1000mm
Powder coated
Tilting writing surface
Pull out keyboard shelf
Storage basket
125mm castors - 2 brake
Option - adjustable height writing surface




AX 841
Modular AV Equipment Trolley
(monitor not included)
680 x 550 x 1500mm
Height to top shelf - 1100mm
Pull out keyboard shelf
100mm castors



AX 840
Modular Equipment Trolley
680 x 550 x 1000mm
Height to top shelf - 1100mm
1 Drawer
100mm castors




AX 636
Theatre Supply Storage Cart
6 Swing out baskets - 2 fixed
Laminated top
4 - 100mm twin castors




This range of equipment is engineered for strength and durability to withstand continuous use year after year. A blend of different materials and finishes are used throughout the range.

All products are fully welded for maximum rigidity and are supported by a full 2 year warranty. A choice of all stainless steel or colour coordinated powder coating is available.



AX 231
Storage Basket and Clamp



AX 232
Vertical Oxy Holder and Clamp




AX 233
IV Handle and Clamp



AX 215
Pogo IV Pole / Stand - Heavy Base
Adjust height by twist mechanism
Height range 1350 - 2100mm
5 Leg base - 75mm castors




Lift lock adjustment



One hand adjustment


AX 2380
Provita Smart IV Stand
One hand operation
Height range 1300 - 2100
Stable five star weighted base
4 hook IV holder
Colour coordinated
Aqua seal / 80mm castors



AX 217
IV Pole / Stand
Height range 1350 - 2050mm
Stainless steel construction
5 Leg weighted base
50mm twin castors



AX 219
IV Pole / Stand - Heavy Base
Height range 1350 - 2050mm
Stainless steel construction
5 Leg weighted base
50mm twin castors



AX 216
IV Pole / Stand Polymer Base
Height range 1350 - 2100mm
5 Leg base



AX 220
IV Transfusion Pump Stand
Height range 1350 - 2100mm
Offset pole for stability
Push handle



AX 528
Oxygen Bottle Trolley G Size
Epoxy coated - 150mm castors
AX 531
Oxygen Bottle Trolley G Size
Stainless steel - 150mm castors



AX 546
Oxygen Bottle Trolley D Size
Epoxy coated - 150mm castors
AX 549
Oxygen Bottle Trolley D Size
Stainless steel - 150mm castors



AX 537
Oxygen Bottle Trolley C Size
Epoxy coated - 150mm castors
AX 540
Oxygen Bottle Trolley C Size
Stainless steel - 150mm castors



AX 990
Charts Trolley - 12 Charts
465 x 530 x 900mm
Stainless steel construction
125mm castors



AX 993
Charts Trolley - 24 Charts
805 x 530 x 900mm
Stainless steel construction
125mm castors
AX 994
Charts Trolley - 36 Charts
1170 x 530 x 900mm



AX 655
Files Trolley
750 x 490 x 900mm
Stainless steel construction
1 Drawer - 125mm deep
Removable clear dividers - Lift up side shelf
125mm castors



AX 524
Multi Stack Gas Cylinder Trolley
Capacity 9 x C size cylinders
Power coated finish



AX 525
Multi Stack Gas Cylinder Trolley
Capacity 12 x C size cylinders
Power coated finish



AX 524SS
Multi Stack Gas Cylinder Trolley
Capacity 9 x C size cylinders
Stainless steel construction



AX 525SS
Multi Stack Gas Cylinder Trolley
Capacity 12 cylinders
6 x C size / 6 x D size
Stainless steel construction



AX 996
Chart Holder A4
To suit Charts Trolleys
310 x 265mm



AX 554
Pan and Bottle Rack - 4
Capacity 4 pans / 5 bottles
Stainless construction
Includes drip tray and back plate



AX 552
Pan and Bottle Rack - 6
Capacity 6 pans / 7 bottles
Stainless construction
Includes drip tray and back plate



AX 501
Bed Cradle
Tubular frame
Stainless steel construction



AX 519
Garbage Trolley-Foot Operated
Foot operated lid action
Stainless construction - PVC lid



AX 177
Lead Apron Rack
700 x 600 x 1420mm
Stainless steel construction
12 apron hooks / corner buffers
125mm castors




AX 175
Garment Rack
1600 x 600 x 1600mm
Stainless steel construction
AX 178 Coat Hangers



AX 534
Gas Bottle Holder
130 x 130 x 350
Hook over / to suit C size
Sheet Metal

AX 543
Gas Bottle Holder
150 dia x 330
Hook over / to suit C size
12mm dia rod to hook into trolley or bed



AX 513
Urinal Bottle Holder
Stainless Steel
AX 510
Urinal Bottle Holder
Nylon Coated



AX 522
Garbage Wall Mount Bracket
310 x 190
Stainless steel
Pre drilled mounting holes

AX 432

Shower Chair

The shower chair offers a three position, drop down armrest, a non-slip retracting footrest fitted with an anti-tilt bracket and comes with a seamless, soft padded seat and backrest. The shower chair's slim line design is made of stainless steel and is fitted with water resistant precision bearing castors which offers excellent manoeuvrability.

Features

- Seat width 440mm
- Seamless padded seat and backrest
- Bed pan holder
- Padded arm rest
- Locking castors
- Stainless steel construction

Options

- Oxygen cylinder holder and clamp
- IV pole and clamp
- Hinged seat
- Safety arms - wrap around

Optional Accessories

(selection shown right)

- AX 436 Hinged Padded Seat
- AX 438 Padded Enclosed Seat
- AX 441 Padded Open Front Seat
- AX 444 Padded Backrest
- AX 462 Oxygen Cylinder Holder & Clamp
- AX 460 IV Pole and Clamp
- AX 455 Safety Arms - Wrap Around
- AX 447 Leg Extension



AX 432SA

Shower Chair

Swing back arms



AX 434

Shower Chair

Seat width 440mm / drop arms
Large rear 150mm castors



AX 433

Shower Chair

Seat width 600mm / drop arms
Weight capacity 175kg



AX 435

Shower Chair - Bariatric

Seat width 660mm / swing back arms
Weight capacity 350kg / heavy duty frame
Reinforced seat / backrest

ACCESSORIES



AX 460

Shower Chair

IV Pole and Clamp



AX 462

Shower Chair

Vertical Oxy Holder and Clamp



AX 455

Shower Chair

Safety Arms - Wrap Around



AX 447

Shower Chair

Leg Extension



AX 708
Fixed Height Patient Trolley
2000 x 685 x 750mm
Stainless steel construction
Removable side rails and end rails
Mediflex mattress
Storage basket / includes oxygen bottle holder
200mm castors / 3 brake / 1 steer



AX 704
Hi Low Mortuary / Concealment Trolley
2200 x 650mm / height range 450 - 1450mm
Stainless steel construction
Hi low function electrically operated by hand controller / battery back up for transportation
Removable side rails and end rails
200mm castors / 3 brake / 1 steer



AX 707
Mortuary / Concealment Trolley Cover

Document Holders
GM document holders offer a variety of configurations for the effective storage of files, binders and documents. The range offers the following options:

- Landscape positioning
- Portrait positioning
- Multiple document holders
- Magnetic mounting option
- Clear or dark tinted options
- Unbreakable polycarbonate construction for some models

Tilt Bins
Modular units with stackable design


- Manufactured using High Impact Polystyrene Plastic (HIPS)
- Transparent bins keep supplies visible
- Supplied with product bin labels for easy identification



AX 9144
Document Holder
Clear - 260 x 254 x 102mm




AX 9142
Document Holder
Clear polycarbonate - 254 x 267 x 51mm




AX 9140
Document Holder
Clear polycarbonate - 445 x 165 x 76mm



AX 9136
Document Holder - Set of 3
Clear polycarbonate - 368 x 165 x 76mm



AX 9146
Document Holder
Clear - 375 x 178 x 102mm




AX 9122
Tilt Bins - 6 Units



AX 9124
Tilt Bins - 5 Units




AX 9126
Tilt Bins - 4 Units



AX 9128
Tilt Bins - 3 Units





The GM range of Maternity products provide for varying conditions of infants, including the ability to bathe. Our range features a twin bassinet which assists in the bonding of the newborns. With Occupational Health and Safety issue being at the forefront we offer 2 types of variable height bassinets units

that offer simple electrically controlled height adjustments with a convenient hand control. This ensures safe, secure handling and treatment of the infant.



AX 609
Single Bassinet Trolley - Adjustable Height
790 x 420mm / height 800 - 1100mm
Electrically operated by hand controller / battery back up
Includes acrylic tub / mattress and 1 x swing out basket



AX 610
Babies Bath - Adjustable Height
790 x 420mm / height 800 - 1100mm
Stainless steel construction
Electrically operated by hand controller / battery back up
Includes acrylic tub with waste outlet connection and hose
1 x swing out basket
AX 611 - Acrylic tub (only) with waste outlet connection and hose




AX 582
Bassinet Trolley - Single
800 x 415 x 910mm
Includes acrylic tub and mattress
2 x white storage trays
Powder coated frame
AX 607 - Stainless steel frame



AX 597
Bassinet Trolley - Twin
800 x 610 x 910mm
Includes acrylic tub and mattress
2 x white storage trays
Powder coated frame
AX 608 - Stainless steel frame



AX 612
Bassinet Trolley - Cabinet
800 x 415 x 910mm
Cabinet with internal adjustable shelf
Includes acrylic tub and mattress



This extensive range of equipment offers handling and transport solutions for linen and other supplies within the healthcare environment. Engineered to withstand heavy loads and frequent use, we have designed these products to be ergonomic

and easy to use. All equipment is fully welded and has a sturdy frame design for maximum durability and is constructed from stainless steel or durable powder coated steel.



AX 300 DOUBLE LINEN SKIP
WITH OPTIONAL HALF HEIGHT SHELF



Foot operated lid







AX 291
Single Linen Skip - foot operated lid
380 x 540 x 850mm
Stainless steel construction
Foot operated PVC lid
75mm castors



AX 300
Double Linen Skip - foot operated lid
780 x 540 x 850mm
Stainless steel construction
Foot operated PVC lid
75mm castors



AX 307
Triple Linen Skip - foot operated lid
1170 x 540 x 850mm
Stainless steel construction
Foot operated PVC lid
75mm castors

Colour Identification Strips
Apply to PVC lid for colour coding (all models)
AX 318R  Red
AX 318B  Blue
AX 318G  Green
AX 318Y  Yellow



AX 319
Soft Close Lid System - Optional
This unique mechanism provides a soft closing action which ensures smooth and quiet operation. This feature can be fitted to all linen skips with PVC lids.



AX 315
Swivel Hooks - Optional
4 hooks per linen bag ring available on all linen skips in the range



AX 303
Triple Linen Skip - no lid
1170 x 540 x 850mm
Stainless steel construction
75mm castors



AX 315
Swivel Hooks - Optional
4 hooks per linen bag ring
available on all linen skips
in the range



AX 294
Double Linen Skip - no lid
780 x 540 x 850mm
Stainless steel construction
75mm castors



AX 286
Single Linen Skip - no lid
380 x 540 x 850mm
Stainless steel construction
75mm castors



AX 320
**Single Linen Skip
Half height shelf
Optional**



AX 733
Moist Laundry Trolley
Including spring loaded base
450 litre capacity



AX 729
Moist Laundry Trolley
Including spring loaded base
320 litre capacity



AX 091
Utility Trolley
720 x 450 x 910mm
Stainless steel frame
2 lift out tubs - 200mm deep
125mm castors



AX 726
Wet and Dry Laundry Trolley
1000 x 650 x 820mm
Including spring loaded base



AX 731
Spring loaded base only



AX 094
Storage Trolley
700 x 450 x 1300mm
Stainless steel frame
4 slide out tubs - 200mm deep
125mm castors

POWDER COATED



AX 307 PC
Triple Linen Skip - foot operated
1170 x 540 x 850mm
Powder coated finish
Foot operated PVC lid
75mm castors



AX 300 PC
Double Linen Skip - foot operated
780 x 540 x 850mm
Powder coated finish
Foot operated PVC lid
75mm castors



AX 291 PC
Single Linen Skip - foot operated
380 x 540 x 850mm
Powder coated finish
Foot operated PVC lid
75mm castors



AX 634
Bulk Collection Trolley
Removable tray
Stainless steel construction
125mm castors



AX 0730
Soiled Linen Trolley
1150 x 575 x 1240mm
Powder coated finish
Gas assisted lid (as shown)
Full bumper surround / 200mm castors





AX 190
Imprest Trolley
1270 x 630 x 1330mm
Powder coated finish
4 Shelves - mesh sides and back
Full bumper surround / 200mm castors
Available with 6 shelves



AX 191
Imprest Trolley
1270 x 630 x 1330mm
Powder coated finish
4 Shelves - enclosed sides and back
Full bumper surround / 200mm castors
Available with 6 shelves



R2530
Residents Clothing Distribution Trolley
Multi compartment trolley designed for distribution of personal linen for residents. Features double sided access. Comprises shelf ledges, vertical rods, hanging rods and push handle.

Sizes available
R2500 - 16 compartment
R2510 - 18 compartment
R2520 - 24 compartment
R2530 - 32 compartment (pictured)



AX 723
Personal Distribution Trolley
960 x 680 x 1310mm
Powder coated frame
20 bays / 20 individual tubs included
Additional sizes available



AX 199
Removable Basket Trolley
750 x 490 x 900mm
Nylon coated baskets - 200mm deep
Push handle / 125mm castors
BASK 0202 Replacement Basket



AX 205
Fixed Basket Trolley
750 x 490 x 900mm
Nylon coated baskets - 200mm deep
Push handle / 125mm castors



AX 147
Clearing Trolley - 2 dished shelves
900 x 490 x 900mm
Stainless steel construction



AX 625
Clearing Trolley
900 x 490 x 900mm
Stainless steel construction
2 break shelves



AX 662
Light Weight Tray Trolley
680 x 550 x 1500mm - 2 x 10 capacity
Tubular stainless steel construction
AX 658 - 1 x 10 capacity
AX 665 - 3 x 10 capacity



AX 150
Clearing Trolley - 3 dished shelves
900 x 490 x 900mm
Stainless steel construction



AX 627
Clearing Trolley
900 x 490 x 900mm
Stainless steel construction
3 break shelves



AX 699
Light Weight Tray Trolley
815 x 480 x 1500mm - 2 x 10 capacity
Depending on tray size
Stainless steel construction



Dished shelf
5mm deep



Break shelf
2mm deep



AX 654
Mobile Preparation Table
1200 x 760 x 900mm
Stainless construction
125mm braking castors

Other size available
AX 651 - 1500 x 760 x 900mm
AX 648 - 1800 x 760 x 900mm
Static Units
AX 645 - 1200 x 760 x 900mm
AX 642 - 1500 x 760 x 900mm
AX 639 - 1800 x 760 x 900mm

Our consulting room range offers an extensive selection of product options to accommodate all types of consulting rooms as well as coordinating with the surrounding décor. The inclusion of fully electrically operated examination couches addresses patient safety for ingress and egress to the

couch and allows flexibility for the physician to treat the patient at a comfortable height.

The options of stainless steel or colour coordinated powder coated equipment is available and accommodates different price levels.

INNOVATION SERIES

The Seers Innovation range of variable height examination couches feature fully electric height adjustment, profiling backrest and footrest and feature extreme low height and generous maximum height.



MG 3470

Innovation 3 Section Couch

- All electric operation
- Height range 450 - 970mm
- Backrest angle 0° / +80°
- Footrest angle +10° / -25°
- Upholstery width 700mm
- Moulded base and acrylic capped frame minimise dirt traps for easy cleaning
- Anti-bacterial powder coated frame
- Anti-microbial, fire retardant upholstery
- Profiled upholstery design
- Clearance beneath frame for mobile hoist
- Large diameter castors (125mm)
- Safe working load 250kg



MG 2470

Innovation 2 Section Couch

- All electric operation
- Height range 450 - 970mm
- Backrest angle -25° / +80°
- Upholstery width 700mm
- Moulded base and acrylic capped frame minimise dirt traps for easy cleaning
- Anti-bacterial powder coated frame
- Anti-microbial, fire retardant upholstery
- Profiled upholstery design
- Clearance beneath frame for mobile hoist
- Large diameter castors (125mm)
- Safe working load 250kg

Optional Accessories



6003
Side Rails



6003P
Side Rail Cushions



6005
Paper Roll Holder

6006
Upper Limb Support



6209
Foot Control 3 Section



6207
Foot Control 2 Section

SX 3008

Sterling 3 Section Couch

- Electric variable height with hand control as standard
- Gas strut assisted adjustable headrest from -25° to +40° positive
- Gas strut assisted adjustable footrest from 0° to +70° positive
- Height range from 450 - 970mm
- Standard width 700mm
- Strong and rigid frame design
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- 100mm castors system with locking brakes
- Safe working load 200kg



SX 2007

Sterling 2 Section Couch

- Electric variable height with hand control as standard
- Gas strut assisted adjustable backrest from -25° to +75° positive
- Height range from 450 - 970mm
- Standard width 700mm
- Strong and rigid frame design
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- 100mm castors system with locking brakes
- Safe working load 200kg





AX 357G
Examination Couch
1810 x 600 x 820mm
Adjustable head section
Tubular frame
EXAM 0354 - Square Tube



AX 385
Examination Couch Heavy Duty
1810 x 760 x 820mm
Adjustable head section
Reinforced frame structure
SWL 250kg

PICTURED WITH OPTIONAL DRAWER UNIT



AX 375G
One Step Stool
Step height 230mm
Epoxy coated - non slip tread




AX 393G
Two Step Stool
Step height 230mm
Epoxy coated - non slip tread



AX 378
One Step Stool with handle rail
Step height 250mm
Epoxy coated - non slip tread
Handle to one end



AX 402RHR - Right
Two Step Stool with hand rail
Step heights 200 / 400mm
Epoxy coated - non slip tread
Hand rail available left or right
AX 402LHR - Left



AX 363
Sliding Arm Board
To suit all models of
Examination Couches



AX 381
Paper Roll Holder
To suit all models of
Examination Couches

**Examination Couch
and Drawer Unit combinations**

- AX 385 DLHS
Left Side Drawer Unit
- AX 385 DRHS
Right Side Drawer Unit

AX 382
3 Drawer Unit (optional)
Modular drawer unit to
suit couch AX 385
Drawers on metal roller
runners

POWDER COATED



AX 321
Instrument Trolley
490 x 490 x 900mm



AX 324
Instrument Trolley - 1 Drawer
490 x 490 x 900mm



AX 327
Instrument Trolley - 2 Drawer
490 x 490 x 900mm



AX 330
Instrument Trolley
750 x 490 x 900mm



AX 333
Instrument Trolley - 1 Drawer
750 x 490 x 900mm



AX 336
Instrument Trolley - 2 Drawer
750 x 490 x 900mm

NOVA



AX 342
Instrument Trolley
490 x 490 x 900mm
Powder coated finish -
Stainless steel shelves



AX 345
Instrument Trolley - 1 Drawer
490 x 490 x 900mm
Powder coated finish
Stainless steel shelves



AX 348
Instrument Trolley - 2 Drawer
490 x 490 x 900mm
Powder coated finish
Stainless steel shelves



AX 349
Instrument Trolley - 3 Drawer
490 x 490 x 900mm
Powder coated finish
Stainless steel shelves

STAINLESS STEEL



AX 033
Instrument Trolley
490 x 490 x 900mm
Stainless steel construction
75mm castors



AX 054
Instrument Trolley - 1 Drawer
490 x 490 x 900mm
Stainless steel construction
75mm castors



AX 066
Instrument Trolley - 2 Drawer
490 x 490 x 900mm
Stainless steel construction
75mm castors



AX 076
Instrument Trolley - 3 Drawer
490 x 490 x 900mm
Stainless steel construction
125mm castors



AX 468
Blood Collection Chair
800 x 1050 x 750mm
Powder coated finish



AX 473
Blood Collection Chair
Stainless steel Sliding Footrest
800 x 1050 x 750mm
Powder coated finish



AX 411
Privacy Screen 3 fold
Panel size 900 x 1400mm
Unfolded 2700 x 1700mm
Powder coated frame
Includes curtains
AX 414CURT
Replacement curtain / 3 fold
AX 423CURT
Replacement curtain / 2 fold



AX 262
Surgeon Stool Hand Operated
Height 530 - 710mm
No backrest



AX 265
Surgeon Stool Hand Operated
Height 530 - 710mm
Backrest



AX 2710
Saddle Stool
Hand operated height adjustment
Height range 530 - 730mm



AX 420
Privacy Screen 2 fold
Panel size 900 x 1400mm
Unfolded 1800 x 1700mm
Powder coated frame
Includes curtains



AX 089
Tub Trolley
490 x 490 x 900mm
4 x Removable tubs
Laminate top
Lift-up side shelf



AX 2382
Examination Light Stand
Multi Purpose
Powder coated finish
5 Star base - 50mm castors




AX 9122
Tilt Bins - 6 Units



AX 9124
Tilt Bins - 5 Units



AX 9126
Tilt Bins - 4 Units




AX 9128
Tilt Bins - 3 Units



AX 9144
Document Holder
260 x 254 x 102mm



AX 9142
Document Holder
254 x 267 x 51mm




AX 9140
Document Holder
445 x 165 x 76mm



AX 9136
Document Holder
368 x 165 x 76mm
Set of 3



AX 9146
Document Holder
375 x 178 x 102mm



Provita Series 1
Clean, modern design, Provita examination lamps have the latest LED technology for economic long term operation.

LED bulbs and the specifically developed reflectors focus clean light where you need it. Optimal heat dissipation utilising a double wall lamp head housing reduces any heat output.

Easy and safe positioning of lamp head with integrated ergonomic handle.




PROVL 100011A
LED Examination Lamp with Gooseneck Arm
Working radius: approx 800mm
Bulb life expectancy: 20,000 hours
Switch at lamp head
Total weight: 2.6kg



PROVL 100031A
LED Examination Lamp with Joint Articulated Arm
Working radius: approx 800mm
Bulb life expectancy: 20,000 hours
Switch at lamp head
Total weight: 2.8kg



PROVL 100081A
LED Examination Lamp with Double-Joint Articulated Arm
Working radius: approx 1100mm
Bulb life expectancy: 20,000 hours
Switch at lamp head
Total weight: 2.8kg



AX 2382
Mobile Examination Lamp with Fixed Height Stand
(lamp not included)

Stable tilt resistant roller base
Base diameter 635mm
Castors 50mm with brake
Steel construction
White powder coat finish

PRODUCT GUIDE



PRODUCT	CODE	PAGE
Airway Trolley	AX 106	13
Anaesthetic Trolley - four drawer	AX 102	05
Anaesthetic Trolley - four drawer	AX 104	05
Anaesthetic Trolley - six drawer	AX 105	05
Arm Table - adjustable height	AX 113	13
Arm Table - fixed height	AX 111	13
AV Equipment Trolley - modular	AX 841	15
Babies Bath - adjustable height	AX 610	25
Basket Trolley - fixed	AX 205	30
Basket Trolley - removable	AX 199	30
Bassinet Trolley - single adjustable height	AX 609	25
Bassinet Trolley - single	AX 582	25
Bassinet Trolley - twin	AX 597	25
Bassinet Trolley - single with cabinet	AX 612	25
Bed Cradle	AX 501	19
Blood Collection Chair	AX 468	36
Blood Collection Chair - sliding footrest	AX 473	36
Bowl Stand - single	AX 211	15
Bowl Stand - single SML	AX 248	15
Bowl Stand - double	AX 214	15
Bulk Collection Trolley	AX 634	29
Charts Trolley - 12 charts	AX 990	19
Charts Trolley - 24 charts	AX 993	19
Charts Trolley - 36 charts	AX 994	19
Chart Holder A4	AX 996	19
Clearing Trolley - 2 break shelves	AX 625	31
Clearing Trolley - 2 dished shelves	AX 147	31
Clearing Trolley - 3 break shelves	AX 627	31
Clearing Trolley - 3 dished shelves	AX 150	31
CSSD Case Cart	AX 692	11
CSSD Kimguard Drying Rack	AX 0681	11
CSSD List Cart	AX 278	11
CSSD Packing Table	AX 0680	11
CSSD Preparation Table - adjustable height	AX 0688	11
CSSD Transportation Trolley	AX 0684	11
CSSD Work Station	AX 0666	10
CSSD Work Station Accessories		
Accessory Shelf	AX 0672	10
CPU Holder	AX 0686	10
Keyboard Shelf	AX 0687	10
Lower Storage Shelf	AX 0689	10
Power Board Outlet 6 plug	AX 0690	10
Sterile Wrap Rack	AX 0685	10
Steri-Peel Holder	AX 0691	10
Distribution Trolley - Clothing	R2530	30
Distribution Trolley - Personal	AX 723	30
Document Holder - 260 x 254 x 102mm	AX 9144	23
Document Holder - 254 x 267 x 51mm	AX 9142	23
Document Holder - 445 x 165 x 76mm	AX 9140	23
Document Holder / set of 3 - 368 x 165 x 76mm	AX 9136	23
Document Holder - 375 x 178 x 102mm	AX 9146	23
Files Trolley	AX 655	19
Emergency Treatment Trolley	AX 115	06
Emergency Treatment Trolley	AX 117	06
Emergency Treatment Trolley	AX 119	06
Equipment Trolley - Modular	AX 840	15
Examination Couch	AX 357G	34
Examination Couch Heavy Duty	AX 385	34
Examination Couch Accessories		
3 Drawer Unit	AX 382	34
Paper Roll Holder	AX 381	34
Sliding Arm Board	AX 363	34
Examination Couch - Electric 2 section	MG 3470	32
Examination Couch - Electric 3 section	MG 2470	32
Examination Couch - Electric 2 section	SX 2007	33
Examination Couch - Electric 3 section	SX 3008	33
Garbag Trolley	AX 519	18
Garbag Wall Mount Bracket	AX 522	19
Garment Rack	AX 175	18
Gas Bottle Trolley - C size single	AX 537	18
Gas Bottle Trolley - C size single stainless	AX 540	18
Gas Bottle Trolley - D size single	AX 546	18
Gas Bottle Trolley - D size single stainless	AX 549	18
Gas Bottle Trolley - G size single	AX 528	18

PRODUCT	CODE	PAGE
Gas Bottle Trolley - G size single stainless	AX 531	18
Gas Bottle Holder - hook over	AX 534	19
Gas Bottle Holder - hook over	AX 543	19
Gas Bottle Trolley - 9 x C size	AX 524	18
Gas Bottle Trolley - 12 x C size	AX 525	18
Gas Bottle Trolley - stainless steel 9 x C size	AX 524SS	18
Gas Bottle Trolley - stainless steel 6 x C / 6 x D size	AX 525SS	18
Imprest Trolley	AX 191	30
Imprest Trolley - mesh	AX 190	30
Instrument Trolley - no rail	AX 001	02
Instrument Trolley - no rail	AX 003	02
Instrument Trolley - no rail	AX 006	02
Instrument Trolley - no rail	AX 009	02
Instrument Trolley - rail 3 sides	AX 033	02
Instrument Trolley - rail 3 sides	AX 036	02
Instrument Trolley - rail 3 sides	AX 039	02
Instrument Trolley - rail 3 sides	AX 042	02
Instrument Trolley - one drawer	AX 054	02
Instrument Trolley - one drawer	AX 057	02
Instrument Trolley - one drawer	AX 060	02
Instrument Trolley - one drawer	AX 063	02
Instrument Trolley - two drawer	AX 066	02
Instrument Trolley - two drawer	AX 069	02
Instrument Trolley - two drawer	AX 072	02
Instrument Trolley - two drawer	AX 075	02
Instrument Trolley - two half drawers	AX 073	03
Instrument Trolley - three drawer	AX 076	03
Instrument Trolley - four drawer	AX 080	03
Instrument Trolley - multi drawer	AX 078	05
Instrument Trolley - Powder coated / stainless drawers		
Instrument Trolley	AX 342	35
Instrument Trolley - one drawer	AX 345	35
Instrument Trolley - two drawer	AX 348	35
Instrument Trolley - three drawer	AX 349	35
Instrument Trolley - Powder Coated		
Instrument Trolley	AX 321	35
Instrument Trolley - one drawer	AX 324	35
Instrument Trolley - two drawer	AX 327	35
Instrument Trolley - one rawer	AX 333	35
Instrument Trolley - two drawer	AX 336	35
ICU Work Station SML	AX 701	15
ICU Work Station LRG	AX 702	15
IV Set Up Trolley	AX 145	06
IV Stand	AX 217	17
IV Stand - heavy base	AX 219	17
IV Stand - polymer base	AX 216	17
IV Stand - POGO heavy base	AX 215	17
IV Stand - transfusion pump	AX 220	17
IV Stand - smart IV	AX 2380	17
IV Stand Accessories		
Handle and Clamp	AX 233	17
Storage Basket and Clamp	AX 231	17
Vertical Oxy Holder and Clamp	AX 232	17
Kick Bucket Stand - combined	AX 240	15
Kick Bucket Stand	AX 241	15
Kick Bucket - 11.4 litre	AX 244	15
Kick Bucket Lid	AX 247	15
Lamp Stand - mobile	AX 2382	37
Lamp Stand - multi purpose	AX 2382	36
Laundry Trolley - Moist 320 litre	AX 729	29
- Spring loaded base only	AX 731	29
Laundry Trolley - moist 450 litre	AX 733	29
Laundry Trolley - wet and dry	AX 726	29
Lead Apron Rack	AX 177	18
LED Examination Lamp - Double-Joint Arm	PROVL100081A	37
LED Examination Lamp - Gooseneck Arm	PROVL100011A	37
LED Examination Lamp - Joint Articulated Arm	PROVL100031A	37
Linen Skip - Double	AX 294	27
Linen Skip - Double Foot Operated Lid	AX 300	27
Linen Skip - Half Height Shelf	AX 320	27
Linen Skip - Single	AX 286	27
Linen Skip - Single Foot Operated Lid	AX 291	27
Linen Skip - Triple	AX 303	27
Linen Skip - Triple Foot Operated Lid	AX 307	27

PRODUCT	CODE	PAGE
Linen Skip Powder Coated		
Linen Skip - single foot operated lid	AX 0291 PC	28
Linen Skip - double foot operated lid	AX 0300 PC	28
Linen Skip - triple foot operated lid	AX 0307 PC	28
Linen Skip Options		
Soft Close Lid System	AX 319	27
Swivel Hooks	AX 315	27
Linen Trolley - soiled	AX 0730	29
Lithotomy Trolley	AX 260	13
Mayo Table - slide base	AX 253	13
Mayo Table - two leg base	AX 252	13
Mayo Table - three leg base	AX 254	13
Mayo Table - four leg base	AX 256	13
Medication Trolley - one drawer	AX 677	07
Medication Trolley - three drawer	AX 676	07
Mortuary / Concealment Trolley	AX 704	22
Mortuary / Concealment Trolley Cover	AX 707	22
Nesting Tables - set of three	AX 021	13
Nesting Table - large	AX 024	13
Nesting Table - medium	AX 027	13
Nesting Table - small	AX 030	13
Pan & Bottle Rack - 4	AX 554	19
Pan & Bottle Rack - 6	AX 552	19
Patient Trolley - fixed height	AX 708	22
Plaster Trolley - two drawer	AX 4990	13
Preparation Table - mobile	AX 654	31
Privacy Screen 2 fold	AX 420	36
Privacy Screen 3 fold	AX 411	36
Privacy Screen replacement curtain / 2 fold	AX 423CURT	36
Privacy Screen replacement curtain / 3 fold	AX 414CURT	36
Shower Chair	AX 432	20
Shower Chair - swing back arms	AX 432SA	20
Shower Chair - fixed rear wheels	AX 434	21
Shower Chair - heavy duty	AX 433	21
Shower Chair - bariatric	AX 435	21
Shower Chair Accessories		
IV Pole and Clamp	AX 460	21
Leg Extension	AX 447	21
Safety Arms - Wrap Around	AX 455	21
Vertical Oxy Holder and Clamp	AX 462	21
Spectra 2 Shelf Trolley	AX S0001	08
Spectra 1 Drawer Trolley	AX S0010	08
Spectra 2 Drawer Trolley	AX S0020	08
Spectra Multi Drawer Trolley	AX S0030	08
Step Stool-one step epoxy coated	AX 375G	34
Step Stool-two step epoxy coated	AX 393G	34
Step Stool-one step epoxy coated/handle rail	AX 378	34
Step Stool-two step epoxy coated/handle rail right	AX 402RHR	34
Step Stool-two step epoxy coated/handle rail left	AX 402LHR	34
Step Stool-one step stainless/aluminium tread	AX 387	14
Step Stool-one step stainless/aluminium tread	AX 390	14
Step Stool-two step stainless/aluminium tread	AX 405	14
Step Stool-one step stainless/alum tread/hand rail	AX 391	14
Storage Trolley	AX 094	29
Storage Unit Mobile - 3 trays	AX M2500	09
Storage Unit Mobile - 5 trays	AX M1400	09
Storage Unit Mobile - 7 trays	AX M8000	09
Storage Unit Mobile - 10 trays	AX M2000	09
Storage Unit Static - 7 trays	AX M8500	09
Surgeons Stool foot operated	AX 267	14
Surgeons Stool foot operated	AX 270	14
Surgeons Stool hand operated	AX 262	14
Surgeons Stool hand operated	AX 265	14
Surgeons Stool - hydraulic	AX 273	14
Surgeon Stool - saddle type	AX 2710	14
Resuscitation Trolley - 4 drawers	AX 107	03
Resuscitation Trolley - standard 5 drawers	AX 108	05
Resuscitation Trolley - 6 drawers	AX 110	05
Theatre Stool	AX 277	14
Theatre Supply Storage Cart	AX 636	15
Tilt Bin - 3 Bin Unit	AX 9128	23
Tilt Bin - 4 Bin Unit	AX 9126	23
Tilt Bin - 5 Bin Unit	AX 9124	23
Tilt Bin - 6 Bin Unit	AX 9122	23

PRODUCT	CODE	PAGE
Tray Trolley	AX 699	31
Tray Trolley - Tubular	AX 662	31
Trays - Sixty Forty System		
50mm deep tray	AX M0601	08
Interlocking divider 600mm	AX M1661	08
Interlocking divider 400mm	AX M1641	08
100mm deep tray	AX M0605	08
Interlocking divider 600mm	AX M1665	08
Interlocking divider 400mm	AX M1645	08
200mm deep tray	AX M0602	08
Interlocking divider 600mm	AX M1662	08
Interlocking divider 400mm	AX M1642	08
Stainless Steel Trolley Accessories		
Drawer insert containers	AX 120	07
Drop down work surface	AX 126	07
IV pole and clamp	AX 129	07
Power board	AX 135	07
Push handle	AX 143	07
Storage panel	AX 132	07
Self closing drawer runner	AX 738	07
Vertical oxy bottle holder and clamp	AX 123	07
Waste bin and bracket	AX 142	07
Tub Trolley	AX 089	36
Urinal Bottle Holder - nylon coated	AX 510	19
Urinal Bottle Holder - stainless steel	AX 513	19
Utility Trolley	AX 091	29

GM Medical Pty Ltd
28-30 Hightech Place
Lilydale VIC 3140
T 03 9735 0111
F 03 9735 2104
gmmedical.com.au

ParagonCare

Paragon Care Group of Companies

