



TABLE OF CONTENTS

Treatment Couches	2
Electric Couches	2
Fixed Height Exam Couches	21
Portable Massage Couches	23
Bedsheets	24
Treatment Couch Accessories & Spare Parts	26
Stainless Steel Trolleys	30
Single Trolleys	30
Double Trolleys	37
No Rail Trolleys	41
Resus & Anaesthesia Trolleys	43
IV Stands	46
Mayo Instrument Trolleys	47
Oxygen Trolleys	48
Linen Skip Trolleys	49
Powder Coated Trolleys	53
Medical Carts	56
Medical Carts	56
Computer Trolleys	65
Lighting	66
Examination Lighting	66
Minor Procedure Lighting	72
X-Ray Viewers	74
Cryotherapy	76
Liquid Nitrogen Flask	76
Liquid Nitrogen Dewars	81
Dewar Accessories	86
Blood Chairs & Stools	89
Blood Collection Chair	89
Round Stools	90
Saddle Stools	91
Step Up Stools	92
Portable Massage Chair	94
Ultrasonic Cleaners	95
Water Distillers	100
Wheelchairs	102
Thermometry	105
Blood Pressure	106
Pulse Oximetry	108
New Products	110

TREATMENT COUCHES

ELECTRIC COUCHES



THREE SECTION MEDICAL TREATMENT COUCH

Our most popular product this three equal section couch with dual gas strut supported sections & smooth hi-lo operation represents great value. Australia's number 1 choice for medical centres for many years.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 5 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

- ET33NB - Navy Blue**
- ET33BL - Black**
- ET33GR - Grey**
- ET33CG - Charcoal Grey**

MAINTENANCE: Please make sure your table is serviced yearly

to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

THREE SECTION BARIATRIC TREATMENT COUCH

Every practice should have at least one heavy duty couch to lift bariatric patients. This three equal section couch with electric backrest and dual gas strut supported footrest boasts a large lifting capacity of 280KG to perform these duties with ease.

Features:

- Electric Hi Lo operation by foot switch
- Electric backrest for effortless operation
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Durable foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head section through electric motor & foot switch from to -40 to +65 Degrees
- Adjustable foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 280KG

ORDER CODES:

ET32NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





TWO SECTION TREATMENT COUCH

One of our most popular tables. This strong two section table will withstand the rigours of patient treatment whether it is used in a demanding Physiotherapy clinic or a small GP practice you can rest assured it will perform. With a sleek design, great features and a foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 67cm
- Leg section length: 123cm
- Adjustable head section through dual gas operated struts from -20 to +75 Degrees
- Height range from 47cm to 92cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 5 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

- ET2NB - Navy Blue**
- ET2BL - Black**
- ET2GR - Grey**

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

THREE SECTION PHYSIO TREATMENT COUCH

The preferred choice of physiotherapists this short head couch provides a strong platform for treatments. With a gas strut supported head and backrest coupled with a smooth hi-lo operation this table represents great value.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 40cm
- Mid section length: 54cm
- Leg section length: 93cm
- Adjustable back section through dual gas operated struts from to +80 Degrees
- Adjustable head section through gas operated struts from -45 to +40 Degrees
- Height range from 47.5cm to 93cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 5 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

ET3NB - Navy Blue

ET3BL - Black

ET3GR - Grey

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





PODIATRY / MULTI PURPOSE CHAIR

A versatile podiatry and multipurpose treatment chair. This strong chair includes many features such as nose hole and plug, flip up armrests & a debris tray to suit the many disciplines of the healthcare and beauty industry.

Features:

- Electric Hi Lo operation by foot switch
- Electric seat tilt operation by foot switch
- Sleek design with rounded frame legs for safety
- Comfortable flip up armrests
- Padded adjustable foot rest
- Head section has face hole and headrest plug
- Debris tray included
- Adjustable feet for uneven surfaces
- 5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 220cm (Flat Table & Including Footrest)
- Width: 91cm (Including Armrests)
- Backrest length: 91cm
- Mid section length: 56cm
- Leg section length: 38cm
- Adjustable backrest through dual gas operated struts to +80 Degrees
- Motorised seat tilt to +25 Degrees
- Height range from 56.5cm to 98cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

PCNB - Navy Blue

PCBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

GYNAECOLOGY TREATMENT COUCH

Everything you will need for women's health and general examinations all in one table. This flexible table includes foot pads, stirrups, debris tray, a fill section & nose hole with plug allowing for a multitude of different examinations & treatments. It features a smooth hi-lo adjustment & electric seat tilt with foot control for optimum manoeuvrability and better patient treatment.

Features:

- Electric Hi Lo operation by foot switch
- Electric seat tilt operation by foot switch
- Sleek design with rounded frame legs for safety
- Footpads & stirrups both included
- Debris tray
- Side rails
- Fill section allowing general use
- Face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 185cm
- Width: 60cm
- Backrest length: 91cm
- Mid section length: 46.5cm
- Leg section length: 42.5cm
- Adjustable backrest through dual gas operated struts to +80 Degrees
- Height range from 64cm to 106cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

GCNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





GYNAECOLOGY PREMIUM TREATMENT COUCH

Our premium gynaecology couch provides for electric positioning of height, tilt and mid backrest through a foot control ensuring practitioner comfort. This couch has been designed with an additional headrest section transforming to a full length couch for most general examinations & minor treatments.

Features:

- Electric Hi Lo operation by foot switch
- Electric seat tilt operation by foot switch
- Electric backrest operation by foot switch
- Debris Tray
- Foot Pads
- Sleek design with rounded frame legs for safety
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 181cm
- Width: 70cm
- Headrest length: 61cm
- Backrest length: 76cm
- Seat section length: 40cm
- Adjustable headrest through dual gas operated struts from -30 to +47 Degrees
- Adjustable motorised backrest from to +71 Degrees
- Adjustable motorised seat tilt up to +18 Degrees
- Height range from 57cm to 118cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

GCPRNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CHIROPRACTIC TABLE

This chiropractic drop table is a versatile cost effective table providing treatment flexibility for your practice. All Cervical, Thoracic, Lumbar and Pelvic drops by foot activation with hand tensioning.



Features:

- Electric Hi Lo operation by foot switch
- Height and width adjustable padded headrest
- Height adjustable & extendable padded footrest
- Sleek design with two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 180-193cm
- Width: 55-59cm
- Headrest length: 40cm
- Tapered sections length: 20.5cm & 22.5cm
- Long section length: 69.5cm
- Footrest Length: 13cm
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable long board to +15 Degrees
- Height adjustable headrest from +25 to -20 degrees
- Height range from 54cm to 87cm in approx 30 seconds
- 5cm wheels
- Padding thickness 8cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

CTNB - Navy Blue

CTBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CHIROPRACTIC PREMIUM TABLE



This Chiropractic table is our new premium table. Boasting a premium column lift eliminating any lifting mechanism. All Cervical, Thoracic, Lumbar and Pelvic drops are operated by foot activation with the tensioning done by hand. This allows easy activation for the practitioner. A premium quality & luxury product.

Features:

- Electric Hi Lo operation by foot switch
- Height and width adjustable padded headrest
- Extendable padded footrest
- Luxurious and comfortable table
- Sleek design with two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 178cm
- Width: 54-65cm
- Headrest height adjustment: 10cm
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable angle of cervical from +25 to -20 Degrees
- Adjustable angle of thoracic to +25 Degrees
- Adjustable angle of lumbar to +30 Degrees
- Adjustable angle of pelvic to +30 Degrees
- Adjustable angle of seat rest to +10 Degrees
- Adjustable angle of backrest to +25 Degrees
- Height range from 60cm to 80cm
- 5cm wheels
- Padding thickness 8cm
- Foam Density: 110g/m²
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

CTPRNB - Navy Blue

CTPRBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CHIROPRACTIC PREMIUM FIXED HEIGHT TABLE

This Chiropractic table is our new premium fixed height table. Very similar to our premium column but is fixed in height. All Cervical, Thoracic, Lumbar and Pelvic drops are operated by foot activation with the tensioning done by hand. This allows easy activation for the practitioner. Overall great value for any clinic.



Features:

- Electric Hi Lo operation by foot switch
- Height and width adjustable padded headrest
- Extendable padded footrest
- Sleek design with two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 173cm
- Width: 62cm
- Headrest length: 40cm
- Adjustable angle of cervical from +15 to -20 Degrees
- Adjustable angle of thoracic to +25 Degrees
- Adjustable angle of lumbar to +30 Degrees
- Adjustable angle of pelvic to +30 Degrees
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Fixed height of 55cm
- 5cm wheels
- Padding thickness 8cm
- Foam Density: 110g/m²
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

CTFIXNB - Navy Blue

CTFIXBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FIVE SECTION TREATMENT COUCH NO POSTURAL DRAINAGE



A popular selection for physiotherapists this five section couch offers strong smooth electric hi-lo operation through footswitch along with a generous range of motion for the head & backrest. This couch will allow correct patient positioning while reducing practitioner strain.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Height adjustable padded armrests
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Headrest length: 43cm
- Mid section length: 51cm
- Backrest length: 96cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from -35 to +70 Degrees
- Height range from 50cm to 96cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

ET51NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FIVE SECTION TREATMENT COUCH WITH POSTURAL DRAINAGE

A popular selection for physiotherapists this five section couch offers strong smooth electric hi-lo and postural drainage operation through footswitch along with a generous range of motion for the head & backrest. This couch will allow correct patient positioning while reducing practitioner strain.



Features:

- Electric Hi Lo operation by foot switch
- Electric postural drainage operation by foot switch
- Sleek design with rounded frame legs for safety
- Height adjustable padded armrests
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Headrest length: 43cm
- Mid section length: 51cm
- Backrest length: 96cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from -35 to +70 Degrees
- Height range from 50cm to 96cm in approx 30 seconds
- Motorised postural drainage adjustment to +25 Degrees
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

ET52NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



NEUROLOGICAL BOBATH TREATMENT COUCH

Designed for rehabilitation this strong two section neurological bobath has been designed with the therapist in mind. The large table top provides plenty of space for patient and therapist and three nose holes and plugs provides flexibility of treatment.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Head section has three face holes and plugs
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 190cm
- Width: 120cm
- Backrest length: 65cm
- Long section length: 121cm
- Adjustable backrest through gas operated strut to +65 Degrees
- Height range from 49cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

NT120NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CONTOURED MASSAGE COUCH

A comfortable and adjustable massage couch with contoured surface to increase patient access. The flexible padded face rest and prone arm rest ensure patient comfort and relaxation during treatment.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Comfortable head cushion with prone arm rest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 195cm
- Width: 70cm
- Headrest length: 32cm
- Long section length: 165cm
- Adjustable headrest through gas operated strut from +30 to -75 Degrees
- Height range from 48cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

CMNB - Navy Blue

CMBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





TILT TABLE TWO SECTION

Our tilt table offers a simple out of the box solution for physiotherapists and in home care patients alike. This two section tilt table including patient straps, activity tray and individually adjustable foot plates will suit the most demanding rehabilitation regime and perform for the long term.

Features:

- Electric hi-lo & tilt operation by hand switch
- Tilt angle inclinometer included
- Activity tray & patient straps included
- Head section has face hole and plug
- Sleek design with rounded frame legs for safety
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 194cm
- Width: 64cm
- Backrest length: 66.5cm
- Long section length: 124cm
- Adjustable backrest through gas operated strut to +75 Degrees
- Adjustable individual foot plates from +20 to -10 Degrees
- Motorised tilt to +85 Degrees
- Height range from 59cm to 104cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

TTNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FOUR SECTION PHYSIO TRACTION TREATMENT COUCH

This traction couch offers a great solution for static and intermittent lumbar and cervical traction treatments. Including a flexion stool and angle adjustable machine mount this couch is an out of the box solution for many practitioners.



Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Upholstered traction stool
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 70cm
- Headrest length: 40cm
- Long seat length: 39cm
- Short seat length: 29cm
- Backrest length: 81cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from +40 to -45 Degrees
- Max traction of 16cm between second and third sections
- Height range from 52cm to 94cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

TRACNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



BEAUTY DAY SPA MASSAGE COUCH

This electric hi lo couch is perfect for use in beauty, day spa and massage. It can be used in numerous comfortable positions with the head and leg section adjusted by gas struts. It has removeable padded arm rests on each side that are also height adjustable. The wings next to the head piece are removeable for patient comfort in massage allowing them to rest their arms on the prone arm rest.

Features:

- Electric Hi Lo operation by foot switch
- Luxurious soft and comfortable head cushion
- Contoured comfortable table top
- Sleek design with rounded frame legs for safety
- Removeable padded wing armrests
- Padded prone arm rest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 186cm
- Width: 50cm
- Backrest length: 58cm
- Mid section length: 45cm
- Legrest length: 41cm
- Long section length: 121cm
- Adjustable backrest through gas operated strut to +85 Degrees
- Adjustable footrest through gas operated strut from +5 to -75 Degrees
- Height range from 57cm to 98cm in approx 30 seconds
- Padding thickness 7.5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

BTW - White

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

THREE SECTION MEDICAL TREATMENT COUCH WITH SIDE RAILS

Our most popular product this three equal section table with side rails & smooth hi-lo operation represents great value.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 5 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

SRET33NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



SINGLE SECTION TREATMENT / CHANGE COUCH WITH SIDE RAILS



This strong single section couch will withstand the rigours of in clinic patient treatment, in a special needs facility or home as a change couch. It has retractable side rails for added patient safety. With a sleek design, great features and a foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot surround bar
- Sleek design with rounded frame legs for safety
- Padded side rails
- Face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 180cm
- Width: 66cm
- Height range from 48cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 5 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007
- Weight Capacity: 225KG

ORDER CODES:

ET1NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FIXED HEIGHT EXAM COUCHES



FREE STANDING TREATMENT COUCH

Our free standing treatment couch is the workhorse of any clinic. Flat packed with a two minute assembly and available in a number of colours choices with and without face hole this couch is a popular choice.

Features:

- Robust
- New lift and lock mechanism for easy 2 minute assembly
- Lift up backrest for assessment
- Face hole and non-face hole models available
- Flat packed to save on freight

Specification:

- Total Length: 195cm
- Width: 71cm
- Height: 82cm
- Weight Capacity: 280KG
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007

ORDER CODES:

- FSBL - BLACK**
- FSBLF - BLACK WITH FACEHOLE**
- FSGR - GREY**
- FSGRF - GREY WITH FACEHOLE**
- FSNB - NAVY BLUE**
- FSNBF - NAVY BLUE WITH FACEHOLE**
- FSNB - WHITE**
- FSNBF - WHITE WITH FACEHOLE**

[Click For More Information](#)

FREE STANDING TREATMENT COUCH WITH DRAWERS

Our three drawer free standing treatment couch is the workhorse of any clinic. Flat packed with a simple & quick assembly this couch is a popular choice.

Features:

Robust
New lift and lock mechanism for easy 2 minute assembly
Lift up backrest for assessment
Non-face hole model only
Three drawer with simple assembly
Flat packed to save on freight

Specification:

Total Length: 195cm
Width: 71cm
Height: 82cm
Drawer: 3
Weight Capacity: 280KG
Padding thickness 5cm
Foam Density: 36kg/m3
Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
Certification: CE, FDA, ISO9001-2008, ISO13485-2007

ORDER CODES:

FSBLDR - BLACK

FSGRDR - GREY

FSNBDR - NAVY BLUE

[Click For More Information](#)



PORTABLE MASSAGE COUCHES



PORTABLE MASSAGE COUCH

Aluminium framed portable massage couch with an auto lock cabling system, rounded corners and a multi density 5cm foam system. A great choice for mobile practitioners.

Features:

- Aluminium frame and legs
- Manufactured without electric welding to join the aluminium
- Auto lock cabling system
- Rounded corners and deluxe no noise handle
- Luxurious PU vinyl
- Adjustable backrest
- Multi density 5cm foam system

Specifications:

- Height Adjustable: 65-80cm
- Length: 185cm
- Width: 66cm
- Folded Size: 93x67x17cm
- Weight: 12KG
- Weight Capacity: 220KG
- Warranty: 2 Years

ORDER CODE:

PMTNB - NAVY BLUE

PMTBL - BLACK

PMCB - CARRY CASE

[Click For More Information](#)

BEDSHEETS

FITTED DISPOSABLE BEDSHEETS BX100

The best form of protection for your examination couch is a disposable fitted bed sheet. Suitable for various treatment couches and stretchers, made to size and packaged in four packs of 25 for easy distribution it is the preferred solution for many healthcare facilities.

Specifications:

2000 x 750mm
Non-woven polypropylene
Available in navy blue & white
40 GSM
Box of 100 (4 x 25)

ORDER CODES:

BSHFITWH - White

BSHFITNB - Navy Blue

[Click For More Information](#)



NON-FITTED DISPOSABLE BEDSHEETS BX100

The best form of protection for your examination couch is a disposable flat bed sheet. Suitable for various treatment couches, made to size and packaged in ten packs of 10 for easy distribution it is the preferred solution for many healthcare facilities.

Specifications:

2400 x 700mm
Non-woven polypropylene
Available in navy blue & white
40 GSM
Box of 100 (4 x 25)

ORDER CODE:

BSHWH - White

BSHNB - Navy Blue

[Click For More Information](#)



TREATMENT COUCH ACCESSORIES & SPARE PARTS



PAPER TOWEL HOLDER FOR ELECTRIC & FIXED HEIGHT COUCHES

A paper towel holder that fixes to electric & fixed height couches.

A paper towel holder that fixes to electric & fixed height couches.

Order Code:

PTH

[Click For More Information](#)

GAS STRUT SHORT FOR ELECTRIC COUCH

Gas strut short for hi-lo couch. Used on many headrest and footrest sections on electric couches.

Gas strut short for hi-lo couch. Used on many headrest and footrest sections on electric couches.

Order Code:

ETGSS

[Click For More Information](#)



GAS STRUT LONG FOR ELECTRIC COUCHES



Gas strut long for electric couches. Used on many backrest sections on electric couches.

Gas strut long for electric couches. Used on many backrest sections on electric couches.

Order Code:

ETGSL

[Click For More Information](#)

HAND CONTROL FOR ELECTRIC COUCH HI-LO

A hand control for activation of the couch Hi-Lo function.

A hand control for activation of the table Hi-Lo function.

Order Code:

ETHC

[Click For More Information](#)



FOOT CONTROL RED DOT FOR ELECTRIC COUCH HI-LO



A foot control used to activate the electric couch hi-lo function.

A foot control used to activate the electric couch hi-lo function.

Order Code:

ETFCRD

[Click For More Information](#)

FOOT CONTROL YELLOW DOT FOR ELECTRIC COUCH TILT

A foot control with yellow dot used to activate the electric couch tilt function. Used on many of our tables such as the Podiatry chair.

A foot control with yellow dot used to activate the electric couch tilt function. Used on many of our tables such as the Podiatry chair.

Order Code:

ETFCYD

[Click For More Information](#)





CONTROL BOX FOR ELECTRIC COUCHES

Control box for electric couches.

Control box for electric couches.

Order Code:

ETCB

[Click For More Information](#)

SPLITTER FOR ELECTRIC COUCHES

A splitter to allow two foot or hand controls to be connected to the table.

A splitter to allow two foot or hand controls to be connected to the table.

Order Code:

ETSP

[Click For More Information](#)



STAINLESS STEEL TROLLEYS

SINGLE TROLLEYS



STAINLESS STEEL TROLLEY NO DRAWER

A high grade stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE: SST00

[Click For More Information](#)

STAINLESS STEEL TROLLEY ONE DRAWER

A high grade single drawer stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

- High grade 304 stainless steel throughout
- One drawer
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:

SST01

[Click For More Information](#)



STAINLESS STEEL TROLLEY TWO DRAWER

A high grade two drawer stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

- High grade 304 stainless steel throughout
- Two drawer
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)



**ORDER CODE:
SST02**[Click For More Information](#)**STAINLESS STEEL TROLLEY THREE DRAWER**

A high grade three drawer stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

- High grade 304 stainless steel throughout
- Three drawer
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

**ORDER CODE:
SST03**[Click For More Information](#)



STAINLESS STEEL TROLLEY FOUR DRAWER

A high grade four drawer stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

- High grade 304 stainless steel throughout
- Four drawer
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE: SST04

[Click For More Information](#)

STAINLESS STEEL TROLLEY ONE DRAWER WITH LOCK

A high grade single stainless steel trolley with one lockable drawer. Has top and bottom side rails and locks on the castor wheels for safety.

Features:

- High grade 304 stainless steel throughout
- One lockable drawer
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- Top locking drawer



95mm lockable castors
Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:
SST01LK

[Click For More Information](#)



STAINLESS STEEL TROLLEY TWO DRAWER WITH TOP LOCKING DRAWER

High grade single stainless steel two drawer trolley. Top drawer with lock. Has top and bottom side rails and locks on the front castor wheels for safety.

Features:

- High grade 304 stainless steel throughout
- Two drawer with top drawer lockable
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- Top locking drawer
- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:
SST02LK

[Click For More Information](#)

STAINLESS STEEL TROLLEY THREE DRAWER WITH TOP LOCKING DRAWER

A high grade single stainless steel three drawers trolley. Top drawer with lock. Has top and bottom side rails and locks on the castor wheels for safety.

Features:

- High grade 304 stainless steel throughout
- Three drawer with top drawer lockable
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- Top locking drawer
- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:

SST03LK

[Click For More Information](#)



STAINLESS STEEL TROLLEY FOUR DRAWER WITH TOP LOCKING DRAWER

A high grade single stainless steel four drawer trolley. Top drawer with lock. Has top and bottom side rails and locks on the castor wheels for safety.

Features:

- High grade 304 stainless steel throughout
- Four drawer with top drawer lockable
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

125mm deep drawer
Top locking drawer
95mm lockable castors
Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:**SST04LK**[Click For More Information](#)

DOUBLE TROLLEYS



STAINLESS STEEL DOUBLE TROLLEY NO DRAWER

A high grade double stainless steel trolley with no drawer. Has top and bottom side rails and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 95mm lockable castors
- Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE: SSTD00

[Click For More Information](#)

STAINLESS STEEL DOUBLE TROLLEY ONE DRAWER (FULL WIDTH)

A high grade double stainless steel trolley with one full width drawer. Has top and bottom side rails and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- One drawer (full width)
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish



Specifications:

125mm deep drawer
95mm lockable castors
Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:
SSTD01

[Click For More Information](#)

**STAINLESS STEEL DOUBLE TROLLEY TWO
DRAWER (FULL WIDTH)**

A high grade double two drawer stainless steel trolley. Has top and bottom side rails and lockable castors for safety.

Features:

High grade 304 stainless steel throughout
Two drawer (full width)
Bottom shelf
Double safety folds
Smooth welding
Soft closing magnetic drawer
German ball bearing drawer runners
Smooth running lockable castors
Top and bottom rails on three sides
Outstanding finish

Specifications:

125mm deep drawer
95mm lockable castors
Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:
SSTD02

[Click For More Information](#)

STAINLESS STEEL DOUBLE TROLLEY TWO DRAWER (SIDE BY SIDE)

A high grade double two drawer stainless steel trolley. Has top and bottom side rails and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Two drawer (side by side)
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:
SSTD22

[Click For More Information](#)



STAINLESS STEEL DOUBLE TROLLEY FOUR DRAWER (TWO OVER TWO)

A high grade double four drawer stainless steel trolley. Has top and bottom side rails and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Four drawer (two over two)
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

125mm deep drawer
95mm lockable castors
Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:
SSTD04

[Click For More Information](#)

NO RAIL TROLLEYS



STAINLESS STEEL TROLLEY NO RAILS 50X50X90

A high grade stainless steel trolley 50cm wide with no draws or rails. Has lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Outstanding finish

Specifications:

- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:
SSTNR50

[Click For More Information](#)

STAINLESS STEEL TROLLEY NO RAILS 60X50X90

A high grade stainless steel trolley 60cm wide with no draws or rails. Has lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Outstanding finish

Specifications:

- 95mm lockable castors
- Dimensions: 60 x 50 x 90cm (LxWxH)

ORDER CODE:



SSTNR60[Click For More Information](#)**STAINLESS STEEL TROLLEY NO RAILS 80X50X90**

A high grade stainless steel trolley 80cm wide with no draws or rails. Has lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Outstanding finish

Specifications:

- 95mm lockable castors
- Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:**SSTNR80**[Click For More Information](#)

RESUS & ANAESTHESIA TROLLEYS



STAINLESS STEEL ANAESTHETIC TROLLEY FOUR DRAWER

A strong high quality anaesthetic trolley with two standard and two deep drawer with rails on the top and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Four drawer
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top rails on three sides
- Outstanding finish

Specifications:

- Drawer Sizes: 2 x 125mm & 2 x 250mm
- 95mm lockable castors
- Dimensions: 50 x 50 x 107cm (LxWxH)

ORDER CODE: SSAT04

[Click For More Information](#)

STAINLESS STEEL ANAESTHETIC TROLLEY SIX DRAWER

A strong high quality six drawer anaesthetic trolley with two side shelves with rails. It has a large top surface with rails on three sides and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Six drawer
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top rails on three sides
- Two side shelves with rails
- Outstanding finish

Specifications:

Drawer Sizes: 1 - 75mm, 2&3 - 150mm, 4 - 75mm, 5&6 - 100mm
 95mm lockable castors
 Dimensions: 90 x 50 x 97cm (LxWxH)

ORDER CODE:
SSAT06

[Click For More Information](#)



STAINLESS STEEL RESUSCITATION TROLLEY FIVE DRAWER

A strong high quality five drawer resuscitation trolley with rails on both shelves and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Five drawer
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top & bottom rails on three sides
- Bottom shelf
- Outstanding finish

Specifications:

Drawer Sizes: 4 x half, 1 x full width
95mm lockable castors
Dimensions: 80 x 50 x 97cm (LxWxH)

ORDER CODE:
SSRT05

[Click For More Information](#)

STAINLESS STEEL RESUSCITATION TROLLEY SEVEN DRAWER

A strong high quality seven drawer resuscitation trolley with rails on top and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Seven drawer
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top & bottom rails on three sides
- Bottom shelf
- Outstanding finish

Specifications:

Drawer Sizes: 6 x 125mm half width, 1 x 250mm full width
95mm lockable castors
Dimensions: 80 x 50 x 97cm (LxWxH)

ORDER CODE:
SSRT07

[Click For More Information](#)



IV STANDS



STAINLESS STEEL IV STAND TWO HOOK

A strong two hook stainless steel IV stand with a 5 leg weighted base and easy mobility.

Specifications:

Two hook
Stainless steel construction
Five Leg weighted base
Smooth rolling and lockable castors
Height range 135-205cm

ORDER CODE:
IVP2H

[Click For More Information](#)

STAINLESS STEEL IV STAND FOUR HOOK

A strong four hook stainless steel IV stand with a five leg weighted base and easy mobility.

Specifications:

Four hook
Stainless steel construction
Five Leg weighted base
Smooth rolling and lockable castors
Height range 135-205cm

ORDER CODE:
IVP4H

[Click For More Information](#)



MAYO INSTRUMENT TROLLEYS



MAYO INSTRUMENT TROLLEY

Strong stainless steel mayo trolley with removeable tray, height adjustment and lockable castors for safety.

Specifications:

304 stainless steel construction
Removable tray: 60 x 40 x 2cm (LxWxD)
Telescopic height adjustment: 85 - 140cm
65mm lockable castors

ORDER CODE:

MIT

[Click For More Information](#)

OXYGEN TROLLEYS



OXYGEN TROLLEY SIZE D

A strong and portable oxygen bottle trolley to suit D size cylinders.

Specifications:

Suits D size cylinder
304 stainless steel throughout
securing chain
rugged wheels
Dimensions: 38 x 27 x 100cm

ORDER CODE: OTD

[Click For More Information](#)

LINEN SKIP TROLLEYS



SINGLE LINEN SKIP

Single stainless steel linen skip with no lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
80mm lockable castors
Dimensions: 40 x 56 x 85cm

ORDER CODE:
SLT

[Click For More Information](#)

DOUBLE LINEN SKIP

Double stainless steel linen skip with no lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
80mm lockable castors
Dimensions: 78 x 56 x 85cm

ORDER CODE:
DLT

[Click For More Information](#)





TRIPLE LINEN SKIP

Triple stainless steel linen skip with no lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
80mm lockable castors
Dimensions: 118 x 56 x 85cm

ORDER CODE:
TLT

[Click For More Information](#)

SINGLE LINEN SKIP WITH LID

Single stainless steel linen skip with foot operated acrylic lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
Foot operated acrylic lid
80mm lockable castors
Dimensions: 40 x 56 x 87cm

ORDER CODE:
SLTL

[Click For More Information](#)





DOUBLE LINEN SKIP WITH LID

Double stainless steel linen skip with foot operated acrylic lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
Foot operated acrylic lid
80mm lockable castors
Dimensions: 78 x 56 x 87cm

ORDER CODE:
DLTL

[Click For More Information](#)

TRIPLE LINEN SKIP WITH LID

Triple stainless steel linen skip with foot operated acrylic lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
Foot operated acrylic lid
80mm lockable castors
Dimensions: 118 x 56 x 87cm

ORDER CODE:
TLTL

[Click For More Information](#)





LINEN SKIP BAG STANDARD

Features:

Constructed from purpose made 300-denier polyester twill fabric known as Durotech or Sharkskin
Heavy Duty Drawcord
Industrial strength cord lock
Bottom handles inside and out
Strong double sewn bottom seams
Four high strength polyester tabs for hanging the bag on a frame
Fits most trolleys
Available in a range of colours

ORDER CODES:

SLBR - Red

SLBB - Blue

SLBW - White

[Click For More Information](#)

LINEN SKIP BAG SOLUBLE SEAM BX250

Features:

Constructed from polyethelene
Dissolvable seam
Detachable water soluble tie
Eliminates handling of soiled or infected laundry
Standard size 710x990mm
Available in a range of colours

ORDER CODES:

SSBR - Red

SSBY - Yellow

SSBP - Purple

[Click For More Information](#)



POWDER COATED TROLLEYS



POWDER COATED TROLLEY NO DRAWER

With smooth castors and an outstanding finish our no drawer powder coated trolley is a great alternative to stainless.

Features:

- Bottom shelf
- Smooth running lockable castors
- Top and bottom side rails on three sides
- Outstanding finish

Specifications:

- 95mm front castors
- Dimensions: 50 x 50 x 96.7cm (LxWxH)
- Powder coated finish

ORDER CODE:
PCT00

[Click For More Information](#)

POWDER COATED TROLLEY ONE DRAWER

With smooth castors and an outstanding finish our single drawer powder coated trolley is a great alternative to stainless.

Features:

- One drawer
- Bottom shelf
- Soft closing magnetic drawer
- Smooth running lockable castors
- Top and bottom side rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm front castors
- Dimensions: 50 x 50 x 96.7cm (LxWxH)
- Powder coated finish



**ORDER CODE:
PCT01**[Click For More Information](#)**POWDER COATED TROLLEY TWO DRAWER**

With smooth castors and an outstanding finish our two drawer powder coated trolley is a great alternative to stainless.

Features:

- Two drawer
- Bottom shelf
- Soft closing magnetic drawer
- Smooth running lockable castors
- Top and bottom side rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm front castors
- Dimensions: 50 x 50 x 96.7cm (LxWxH)
- Powder coated finish

**ORDER CODE:
PCT02**[Click For More Information](#)**POWDER COATED TROLLEY THREE DRAWER**

With smooth castors, an outstanding finish our three drawer powder coated trolley is a great alternative to stainless.

Features:

- Three drawer
- Bottom shelf
- Soft closing magnetic drawer
- Smooth running lockable castors
- Top and bottom side rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer



95mm front castors
Dimensions: 50 x 50 x 96.7cm (LxWxH)
Powder coated finish

ORDER CODE:
PCT03

[Click For More Information](#)



POWDER COATED TROLLEY FOUR DRAWER

With smooth castors and an outstanding finish our four drawer powder coated trolley is a great alternative to stainless.

Features:

- Four drawer
- Soft closing magnetic drawer
- Smooth running lockable castors
- Top and bottom side rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm front castors
- Dimensions: 50 x 50 x 96.7cm (LxWxH)
- Powder coated finish

ORDER CODE:
PCT04

[Click For More Information](#)

MEDICAL CARTS



EMERGENCY TROLLEY

This rugged portable crash cart is constructed with stainless steel & an ABS modular surface. With a heavy load rating and complete with accessories it's a must for many healthcare facilities.

Note: Product pictured with central & seals lock. It can only come with one or the other.

Features:

- Five drawer
- Extension table
- CPR board
- Defibrillator Shelf
- Stainless steel IV pole
- Waste bin
- 304 stainless steel wire basket
- Oxygen tank holder
- External power holder
- ABS drawer dividers
- Central or seals lock
- Noise resistant castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Smooth welding

Specifications:

- Dimensions: 650 x 480 x 980mm
- Drawer Dimensions: 550 x 400mm
- Drawer Depths: 1-3 110mm, 4 - 135mm, 5 - 200mm
- Drawer load capacity: 50kg
- Overall load capacity: 200kg
- Weight: 55kg

All accessories included

ORDER CODE:
EMT

[Click For More Information](#)

MEDICINE TROLLEY

For the safe transport of medications and supplies our medicine trolley is the perfect fit. Including all accessories and boasting a number of features this product represents great value.

Features:

- Four drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

- Dimensions: 650 x 480 x 900mm
- Drawer Dimensions: 550 x 400mm
- Drawer Depths: 1-2 120mm, 3 - 140mm, 4 - 200mm
- Top loading: 75kg
- Weight: 42kg

All accessories included

ORDER CODE:
MEDT

[Click For More Information](#)





TREATMENT TROLLEY ONE DRAWER

This single drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

Features:

- One drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

- Dimensions: 650 x 480 x 900mm
- Drawer Depth: 110mm
- Top loading: 75kg
- Weight: 35kg

All accessories included

ORDER CODE:
TTR01

[Click For More Information](#)

TREATMENT TROLLEY TWO DRAWER

This two drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

Features:

- Two drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:



Dimensions: 650 x 480 x 900mm
Drawer Depth: 110mm
Top loading: 75kg
Weight: 35kg

All accessories included

**ORDER CODE:
TTR02**

[Click For More Information](#)



TREATMENT TROLLEY THREE DRAWER

This three drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

Features:

- Three drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

Dimensions: 650 x 480 x 900mm
Drawer Depth: 110mm
Top loading: 75kg
Weight: 35kg

All accessories included

**ORDER CODE:
TTR03**

[Click For More Information](#)

ANAESTHESIA TROLLEY

Anesthesia trolley with an extension table, waste bin, stainless steel wire basket and soft closing drawers.

Features:

- Six drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Storage rack

Specifications:

Dimensions: 650 x 480 x 980mm
 Drawer Depth: 4 x 110mm, 1 x 140mm
 Weight: 55kg

All accessories included

ORDER CODE:
ANT

[Click For More Information](#)



DRESSING CHANGE TROLLEY

A high quality two drawer trolley for the changing of dressings.

Features:

- Two drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Two mounted buckets

Specifications:

Dimensions: 650 x 480 x 900mm
Drawer Depth: 110mm
Weight: 35kg

All accessories included

ORDER CODE:
DCTR

[Click For More Information](#)

INFUSION TROLLEY

A fully featured infusion trolley including bins, basket & adjustable height IV bag pole.

Features:

- One drawer
- Adjustable height rack
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

Dimensions: 650 x 480 x 900mm
Drawer Depth: 110mm
Weight: 22kg

All accessories included

ORDER CODE:
INFTR

[Click For More Information](#)





INSTRUMENT TROLLEY NO SHELF

A strong ABS instrument trolley for transport of goods through healthcare facilities.

Features:

- Luxurious push handle
- Strong polymer shelves
- Easy to clean
- Castors with brakes

Specifications:

Dimensions: 650 x 480 x 900mm
Weight: 22kg

ORDER CODE:
INST00

[Click For More Information](#)

INSTRUMENT TROLLEY ONE SHELF

A strong ABS instrument trolley with a single shelf for transport of goods through healthcare facilities.

Features:

- Luxurious push handle
- Strong polymer shelves
- Easy to clean
- Castors with brakes

Specifications:

Dimensions: 650 x 480 x 900mm
Weight: 22kg

ORDER CODE:
INST01

[Click For More Information](#)





DIFFICULT AIRWAY TROLLEY

Six drawer difficult airway trolley with extension table, luxurious push handle and all accessories included.

Features:

- Six drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Lockable end cabinet

Specifications:

Dimensions: 650 x 480 x 1070mm
 Drawer Depth: 3 x 110mm, 2 x 135mm, 1 x 195mm
 Weight: 60kg

All accessories included

ORDER CODE:
DATR

[Click For More Information](#)

MEDICAL RECORDS TROLLEY NO DRAWER

A full featured medical records trolley including all accessories. A great addition to any healthcare facility.

Features:

- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- CT file container

Specifications:

Dimensions: 650 x 480 x 900mm
 Weight: 35kg



All accessories included

**ORDER CODE:
MEDRT0**

[Click For More Information](#)



MEDICAL RECORDS TROLLEY ONE DRAWER

A full featured single drawer medical records trolley including all accessories. A great addition to any healthcare facility.

Features:

- One drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- CT file container

Specifications:

- Dimensions: 650 x 480 x 900mm
- Drawer Depth: 110mm
- Weight: 35kg

All accessories included

**ORDER CODE:
MEDRT1**

[Click For More Information](#)

COMPUTER TROLLEYS



MEDICAL COMPUTER TROLLEY WITH PRINTER SHELF

A great addition to any healthcare facility or hospital. Our versatile workstation trolley is packed full of features and comes with all accessories included.

Features:

- Display arm
- Disinfectant Holder
- Laptop box with lock & key (separate keys for each side)
- Laptop sits inside laptop box underneath worktop and connects to external screen
- Gas Lift Handle for height adjustment
- Slide keyboard holder
- Pull-out mouse tray
- Mouse holder
- Printer tray
- Twin wheel braked castors
- Manual Gas Spring height adjustment

Specifications:

- Length: 54cm
- Width: 70cm
- Height range: 119-142cm
- Work surface height: 81-116cm
- Work surface size: 53x48cm
- Printer Tray Size: 51x68cm
- Castors: 100mm twin wheel braked

Order Code:

MWT

[Click For More Information](#)

LIGHTING

EXAMINATION LIGHTING



PML1 LED EXAMINATION LIGHT

The PML1 LED examination light is a great all-round diagnostic light. It provides a soft white light for better colour rendition and easier diagnosis and treatment. With a small light head and convenient handle with integrated on/off switch it is easy to move without blocking the vision of the operator. The flexible gooseneck allows positioning in a diverse range of angles for optimal illumination of the work area. A stepless dimmer control allows the operator to set the precise illumination required for each use without being constrained to factory preset levels and the heat free LED provides bright white light for up to 50,000 hours of operation.

It is designed for use in Consulting Rooms, Recovery Areas, Emergency Rooms, Gynaecology Examinations and Minor Surgical Procedures.

Features:

- Powerful long-life LED with up to 50,000 hours of operation
- Low energy consumption
- IR and heat free cold white light
- Comfortable guiding handle with integrated On/Off switch
- Cooling through a heatsink so no fans or moving parts
- Superb colour rendering
- Mobile base and wall mount provided with every light

Specifications:

- Light Source: 3W LED
- Light Intensity: 89,000 Lux @ 30cm / 40,000 Lux @ 50cm
- Light Field Diameter: 80-150mm
- Colour Temperature: 4,500-5,000K
- Lamp Head Diameter: 75mm
- Lamp Head Length: 110mm
- Base Diameter: 530mm
- Base Height: 790mm
- Gooseneck Length: 760mm
- Total Height: 1,830mm

Order Codes:

- PML1WH - WHITE**
- PML1BL - BLACK**

[Click For More Information](#)

PML2 LED EXAMINATION LIGHT

The PML2 LED examination light has 6 powerful LED's offering powerful illumination and great value. It provides white light for better colour rendition and easier diagnosis and treatment. With a small light head with integrated handle it is easy to move without blocking the vision of the operator. The flexible goose neck allows positioning in a diverse range of angles for optimal illumination of the work area. A stepless dimmer control allows the operator to set the precise illumination required for each use without being constrained to factory preset levels and the heat free LED provides bright white light for up to 50,000 hours of operation.

It is designed for use in Consulting Rooms, Recovery Areas, Emergency Rooms, Gynaecology Examinations & Minor Surgical Procedures.

Features:

- Powerful long-life LED with up to 50,000 hours of operation
- Low energy consumption
- IR and heat free cold white light
- Small light head with integrated handle
- Cooling through a heatsink so no fans or moving parts
- Superb colour rendering
- Mobile base and wall mount provided with every light


Specifications:

- Light Source: 6W LED
- Light Intensity: 78,000 Lux @ 30cm / 40,000 Lux @ 50cm / 20,000 Lux @ 100cm
- Light Field Diameter: 70-200mm
- Colour Temperature: 4,500-5,500K
- Lamp Head Diameter: 80mm
- Base Diameter: 530mm
- Base Height: 790mm
- Gooseneck: Length 760mm
- Total Height: 1,815mm

Order Codes:

- PML2WH - WHITE**
- PML2BL - BLACK**

[Click For More Information](#)



FLEXLED EXAMINATION LIGHT

FlexLED is one of the most economical and flexible LED examination lights ever made in Australia. Perfect for general practices and workstations.

FlexLED offers over 40,000 lux @ 1 metre, a 4000K colour temperature and unlimited range of movement with a hygienically sealed flexible arm and head. It's so flexible it allows accurate positioning of the lamp head above almost any work surface.

Specifications:

- 40,000+ lux @ 1 metre
- colour temperature 4,000K
- colour rendering (Ra): 95
- light beam angle: 10°
- hygienically sealed arm and head
- unlimited range of movement
- total wattage including driver 12W
- desk, wall and trolley mountable

Applications:

General examination
Nursing and work stations
First aid stations

Order Codes:

Mobile: LEDFLEX + LEDFLEXTM

Wall: LEDFLEX + LEDFLEXWM

Permanent Desk Mount: LEDFLEX + LEDFLEXDMP

[Click For More Information](#)

ULED EXAMINATION LIGHT

ULED's incredible design story is steeped in 60 years medical lighting manufacture experience combined with a decade of LED lighting innovation. ULED brings examination lighting into the future. Cost-effective yet built to last, ULED's ingenious floating spring arm enables intuitive and perfect positioning. Its powerful yet cool beam provides 40,000+ lux @1 metre.

ULED makes cleaning and operation easy with a fully sealed and touch sensitive membrane switch complete with 4 stage dimming. Plus its push-to-release detachable light head makes upgrade or maintenance possible without a technician. ULED is effortless light – just as it should be.

Designed and manufactured in Australia it comes with an impressive 5 year warranty with industry-leading support, and is available in a range of mount options.

Specification:

- 40,000+ lux @ 1 metre
- quick release head
- IP65 light head shaped for easy cleaning/sterilisation
- total power use 14W
- low operational cost
- next-to-zero maintenance
- all metal construction for durability
- mounting brackets for every application (ceiling, wall, extended wall and trolley mount)
- large task area
- generous 2200mm range of movement
- minimum intrusion into the workspace
- retrofittable to UDP Planet Halogen Light mounts
- light weight – less than 5kg including ceiling mounts.
- designed and manufactured in Australia
- complies with AS/NZ IEC60601-2-41
- 5 year warranty

Applications:

Examination
Over Bed
Intensive care
Recovery rooms
Dermatology
Aged Care
First Aid Rooms

Order Codes:

ULEDUC - Ceiling Mount 2.7M

ULEDUCEPS - Ceiling Mount 2.7-3M

ULEDUCECM - Ceiling Mount 3.6M+

ULEDUW - Wall Mount

[Click For More Information](#)





ULED PLUS EXAMINATION LIGHT

ULED's incredible design story is steeped in 60 years medical lighting manufacture experience combined with a decade of LED lighting innovation. ULED Plus brings examination lighting into the future. Cost-effective yet built to last, ULED Plus ingenious floating spring arm enables intuitive and perfect positioning. Its powerful yet cool beam provides 50,000+ lux @1 metre.

ULED Plus makes cleaning and operation easy with a fully sealed and touch sensitive membrane switch complete with 4 stage dimming. Plus its push-to-release detachable light head makes upgrade or maintenance possible without a technician. ULED Plus is effortless light – just as it should be.

Designed and manufactured in Australia it comes with an impressive 5 year warranty with industry-leading support, and is available in a range of mount options.

Specification:

- 50,000+ lux @ 1 metre
- quick release head
- IP65 light head shaped for easy cleaning/sterilisation
- total power use 20W
- low operational cost
- next-to-zero maintenance
- all metal construction for durability
- mounting brackets for every application (ceiling, wall, extended wall and trolley mount)
- large task area
- generous 2200mm range of movement
- minimum intrusion into the workspace
- retrofittable to UDP Planet Halogen Light mounts
- light weight – less than 5kg including ceiling mounts.
- designed and manufactured in Australia
- complies with AS/NZ IEC60601-2-41
- 5 year warranty

Applications:

Examination
Over Bed
Intensive care
Recovery rooms
Dermatology
Aged Care
First Aid Rooms

Order Codes:

ULEDUCPLUS - Ceiling Mount 2.7M

ULEDUWPLUS - Wall Mount

[Click For More Information](#)

MINOR PROCEDURE LIGHTING



PHANTOM LED MINOR PROCEDURE LIGHT

Delivering a cool and efficient 70,000+ lux @ 1 metre, Phantom's sleek, lightweight design and large task area action zone make it perfect for a range of applications – and with only minimal intrusion to the task area. Plus, Phantom's completely sealed IP65 head gives quick access to ON/OFF and 4 stage dimming control, providing easy operation. Plus, its quick-release head is fully removable for technician-free maintenance.

Phantom offers an improved working environment, significantly reducing on-going maintenance with low operational costs and is a "greener" alternative to similar products. Its easy clean surfaces and removable handle allow superior infection control and its precise arm with multiple rotation points give you light positioning freedom and comes complete with an adjustable floating balance.

Designed to last, Phantom is 100% designed and manufactured in Australia for your confidence. The Planet Medical Phantom is the solid state medical light that needs to be experienced.

Specifications:

- powerful 70,000+ lux @ 1 metre
- quick release head
- IP65 light head shaped for easy cleaning/sterilisation
- total power use 31W
- low operational cost and minimal maintenance
- all metal construction for durability
- large task area
- extensive 3650mm range of movement
- minimum intrusion into the workspace
- multi axis head rotation
- 3 continuous rotation points
- light weight – less than 11kg including ceiling mounts.
- designed and manufactured in Australia
- complies with AS/NZ IEC60601-2-41 (Medical electrical equipment – Particular requirements for safety – Surgical luminaires and luminaires for diagnosis)
- prompt supply of spare parts, technical service and support
- 5 year warranty

Applications:

Examination
Minor Procedure
Over Bed
Intensive Care
Recovery Rooms
Dermatology
Aged Care
First Aid Rooms

Order Codes:

LPHANC: Ceiling Mount 2.7M

LPHANW: Wall Mount

[Click For More Information](#)

X-RAY VIEWERS



SLIMLINE X-RAY VIEWERS

A high quality slimline LED xray viewer with 2.4cm thickness, mitsubishi screen, high luminance, 100000 hours working time and very simple operation. A great addition to any clinic.

Features:

- Mitsubishi Screen
- 2.4cm Thickness
- No Globes Required
- Lightweight
- High Luminance

Specifications:

- 10000K Colour Temperature
- 100,000 Hours Working Time
- 4000cd/m2 / Average 13000 Lux Luminance

Dimensions (L x W x H): 47x2.4x54.5cm (Single), 85x2.4x54.5cm (Double), 123x4.1x54.5cm (Triple), 161x4.1x54.5cm (Quad)

Order Codes:

- SBXSL - Single Bay**
- DBXSL - Double Bay**
- TBXSL - Triple Bay**
- QBXSL - Quad Bay**

[Click For More Information](#)

STANDARD X-RAY VIEWERS

A high quality lightweight LED xray viewer with a mitsubishi screen, high luminance and 100000 hours working time. Great for any clinic.

Features:

Mitsubishi Screen
12cm Thickness
LED Illumination
Lightweight

Specifications:

10000K Colour Temperature
100,000 Hours Working Time
2000cd/m2 / Average 6300 Lux Luminance

Dimensions: 44x12x59cm (Single), 80x12x59cm (Double),
116x12x59cm (Triple), 150x12x59cm (Quad)

Order Codes:

SBXST - Single Bay

DBXST - Double Bay

TBXST - Triple Bay

QBXST - Quad Bay

[Click For More Information](#)



CRYOTHERAPY

LIQUID NITROGEN FLASK



NITROSPRAY CRYOTHERAPY FLASK

Premier Nitrospray cryotherapy flask constructed of a double-walled vacuum stainless steel container to hold liquid nitrogen. The outflow track extension port, located on the top of the instrument has a luer-lock connector to easily engage & disengage spray tips and other attachments as required.

Features:

- Double walled vacuum cryotherapy flask
- Easy connect luer lock connector
- Relief valve precisely regulates pressure
- Curved trigger provides unobstructed view of treatment site
- Easy to hold grooved head
- Durable stainless steel construction
- Spray tip protector prevents accidental damage and contact with patient

Specification:

- Holding Time: 6-8 Hours
- Intermittent Use: 3-6 Hours Treatment Time
- Package Inclusions: Flask, Set of five spray tips & Spray tip protector
- Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)
- Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

ORDER CODES:

LNSLITEPLUS - 300ml

LNSPLUS - 475ml

[Click For More Information](#)

CRYOTHERAPY PACKAGE ONE

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle, Extractor & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle, Extractor & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal



Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

Order Code:

LNPACK1

[Click For More Information](#)

CRYOTHERAPY PACKAGE TWO

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal



Dewar Specification:

Capacity: 20L
Weight (empty): 9.5kg
Caliber (neck opening of dewar): 50mm
Outer diameter of dewar: 394mm
Height of dewar: 672mm
Static evaporation rate: 0.12L/d
Static holding time: 168 days
Ladle outer diameter: 38mm
Ladle height: 120/276mm
No. ladles per dewar: 3~6

Order Code:**LNPACK2**[Click For More Information](#)**NITROSPRAY CRYOTHERAPY SPRAY TIPS**

Spray tips suitable for the Nitrospray cryotherapy flasks.

Features:

Suitable for use with the Nitrospray 300ml & 475ml flasks
Easy connect luer lock connector
Comes in 16, 17, 18, 19 & 20 gauge sizes
Available individually or as a set of five

ORDER CODES:**LNST16G - 16 Gauge****LNST17G - 17 Gauge****LNST18G - 18 Gauge****LNST19G - 19 Gauge****LNST20G - 20 Gauge****LNSTSET - 16, 17, 18, 19 & 20 Gauge**[Click For More Information](#)

NITROSPRAY SPRAY TIP PROTECTOR



A protector for the spray tip of the Nitrospray cryotherapy flask.

A protector for the spray tip of the Nitrospray cryotherapy flask.

Order Code:

LNSTP

[Click For More Information](#)

NITROSPRAY ACCESSORY ADAPTOR FOR SCREW TIPS

An adaptor that allows the connection of screw tips from brands such as Brymill.

An adaptor that has a luer lock connection on one end and a screw tip on the other. This allows the connection of spray tips from brands such as Brymill.



Order Code:

LNSABA

[Click For More Information](#)

NITROSPRAY SPRAY TIP EXTENSION BENT



A long curved needle spray tip which connects to the luer lock connection.

A long curved needle spray tip which connects to the luer lock connection. Allows better patient access when treating difficult areas.

Order Codes:

LNSTEB

[Click For More Information](#)

NITROSPRAY NEEDLE EXTENDER RIGID

A rigid extender for the Nitrospray spray tips.

A rigid steel extender for the existing luer lock connection. Allows the tip better reach and patient access when treating difficult areas.



Order Code:

LNSNER

[Click For More Information](#)

LIQUID NITROGEN DEWAR



10L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 10L dewar with a safety lock lid and low liquid nitrogen consumption.

Features:

- Made of a high-strength aluminium alloy
- Advanced powder coating technology
- Designed for low liquid nitrogen consumption
- 5 year vacuum warranty
- Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

Specifications:

- Capacity: 10L
- Weight (empty): 6.0kg
- Caliber (neck opening of dewar): 50mm
- Outer diameter of dewar: 300mm
- Height of dewar: 552mm
- Static evaporation rate: 0.12L/d
- Static holding time: 86 days
- Ladle outer diameter: 38mm
- Ladle height: 120mm
- No. ladles per dewar: 3~6
- Warranty: 5 Years Vacuum, 2 Years Container

Straws holding capacity:
Single 120mm (0.5mL): 792pc
Single 120mm (0.25mL): 1788pc

ORDER CODE: LNC10

[Click For More Information](#)

20L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 20L dewar with a safety lock lid and low liquid nitrogen consumption.

Features:

Made of a high-strength aluminium alloy
Advanced powder coating technology
Designed for low liquid nitrogen consumption
5 year vacuum warranty
Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

Specifications:

Capacity: 20L
Weight (empty): 9.5kg
Caliber (neck opening of dewar): 50mm
Outer diameter of dewar: 394mm
Height of dewar: 672mm
Static evaporation rate: 0.12L/d
Static holding time: 168 days
Ladle outer diameter: 38mm
Ladle height: 120/276mm
No. ladles per dewar: 3~6

Straws holding capacity:
Single 120mm (0.5mL): 792pc
Single 120mm (0.25mL): 1788pc
Double 276mm (0.5mL): 1284
Double 276mm (0.25mL): 2832

ORDER CODE: LNC20

[Click For More Information](#)





30L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 30L dewar with a safety lock lid and low liquid nitrogen consumption.

Features:

- Made of a high-strength aluminium alloy
- Advanced powder coating technology
- Designed for low liquid nitrogen consumption
- 5 year vacuum warranty
- Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

Specifications:

- Capacity: 30L
- Weight (empty): 12.9kg
- Caliber (neck opening of dewar): 50mm
- Outer diameter of dewar: 462mm
- Height of dewar: 706mm
- Static evaporation rate: 0.12L/d
- Static holding time: 254 days
- Ladle outer diameter: 38mm
- Ladle height: 120/276mm
- No. ladles per dewar: 3~6

- Straws holding capacity:
 - Single 120mm (0.5mL): 792pc
 - Single 120mm (0.25mL): 1788pc
 - Double 276mm (0.5mL): 1284
 - Double 276mm (0.25mL): 2832

ORDER CODE: LNC30

[Click For More Information](#)

CRYOTHERAPY PACKAGE ONE

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle, Extractor & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle, Extractor & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal



Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

Order Code:

LNPACK1

[Click For More Information](#)

CRYOTHERAPY PACKAGE TWO

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal



Dewar Specification:

Capacity: 20L
Weight (empty): 9.5kg
Caliber (neck opening of dewar): 50mm
Outer diameter of dewar: 394mm
Height of dewar: 672mm
Static evaporation rate: 0.12L/d
Static holding time: 168 days
Ladle outer diameter: 38mm
Ladle height: 120/276mm
No. ladles per dewar: 3~6

Order Code:**LNPACK2**[Click For More Information](#)

DEWAR ACCESSORIES



EXTRACTOR FOR 10L DEWAR

A safe way for extraction of liquid nitrogen from a 10L dewar.

Extractor to suit 10 litre dewar.

Order Code:

LNEXT10

[Click For More Information](#)

EXTRACTOR FOR 20L & 30L DEWARS

A safe way for extraction of liquid nitrogen from a 20 & 30L dewar.

Extractor to suit 20 and 30 litre dewars.

Order Code:

LNEXT2030

[Click For More Information](#)





LADLE FOR 10L DEWARS

A ladle to remove liquid nitrogen from a 10L dewar.

Liquid nitrogen ladle to suit 10 litre dewars.

Order Code:

LNL10

[Click For More Information](#)

LADLE FOR 20L & 30L DEWARS

A ladle to remove liquid nitrogen from a 20 & 30L dewar.

Liquid nitrogen ladle to suit 20 and 30 litre dewars.

Order Code:

LNL20

[Click For More Information](#)





LIQUID NITROGEN GLOVES

A necessity for every practice handling liquid nitrogen is a high protection glove.

Protective gloves for use when decanting liquid nitrogen from a dewar to flask.

Order Code:

LNGLV

[Click For More Information](#)

TROLLEY FOR 20L & 30L DEWARS

A trolley that allows the safe transport of a liquid nitrogen dewar through healthcare facilities.

Rolling trolley to suit 20 & 30L dewars.

Order Code:

LNTR

[Click For More Information](#)



BLOOD CHAIRS & STOOLS

BLOOD COLLECTION CHAIR



BLOOD COLLECTION CHAIR

Comfortable and supportive blood collection chair including arm rests and a foldable foot rest.

Specifications:

Comfortable soft, supportive high density foam
Extra wide Seat: 55cm
Width including arm rests: 91.5cm
Arm rests width: 18.5cm
Length: 65cm
Length including foot rests: 84cm
Foldable foot rest
Weight capacity: 175kg

ORDER CODE:

BCCGR - GREY

BCCBL - BLACK

[Click For More Information](#)

ROUND STOOLS



ROUND STOOL STANDARD

A great value round stool suitable for many applications and practices.

Specifications:

Adjustable height: 49cm - 63cm
Cushion diameter: 35cm
Foam thickness: 6cm

ORDER CODES:

RSNB - Navy Blue

RSBL - Black

RSGR - Grey

[Click For More Information](#)

ROUND STOOL PREMIUM WITH FOOT RING

A premium round stool with a height adjustable foot ring and strong base.

Specifications:

Adjustable height: 49cm - 63cm
Cushion diameter: 35cm
Foam thickness: 6cm
Base: 56cm
Adjustable height foot ring

ORDER CODES:

PRSNB - Navy Blue

PRSBL - Black

PRSGR - Grey

[Click For More Information](#)



SADDLE STOOLS



SADDLE STOOL STANDARD

A great value saddle stool suitable for many applications and practices.

Specifications:

Adjustable height: 49cm - 63cm

Cushion diameter: 35cm

Foam thickness: 6cm

ORDER CODES:

SSNB - Navy Blue

SSBL - Black

SSGR - Grey

[Click For More Information](#)

SADDLE STOOL PREMIUM WITH FOOT RING

A premium saddle stool with a height adjustable foot ring and strong base.

Specifications:

Adjustable height: 49cm - 63cm

Cushion diameter: 35cm

Foam thickness: 6cm

Base: 56cm

Adjustable height foot ring

ORDER CODES:

PSSNB - Navy Blue

PSSBL - Black

PSSGR - Grey

[Click For More Information](#)



STEP UP STOOLS

SINGLE STEP STOOL



Perfect for any examination couch this sturdy single step stool is made from 1" stainless steel tubing with reinforced rubber tips & a non slip ribbed rubber platform.

Specifications:

Stainless steel frame
Non-slip ribbed rubber platform
Reinforced rubber tips
Manufactured with 1" steel tubing
Dimensions: 455x245x200 (LxWxH)
Weight Capacity: 150KG

ORDER CODE: SST

[Click For More Information](#)

DOUBLE STEP STOOL

Perfect for any examination couch this sturdy double step stool is made from 1" stainless steel tubing with reinforced rubber tips & a non slip ribbed rubber platform.

Specifications:

Stainless steel frame
Non-slip ribbed rubber platform
Reinforced rubber tips
Manufactured with 1" steel tubing
Overall Dimensions: 400x550x470mm (LxWxH)
First Step: 250mm
Second Step: 220mm
Weight Capacity: 150KG

ORDER CODE: DST

[Click For More Information](#)



PORTABLE MASSAGE CHAIR



PORTABLE MASSAGE CHAIR

Perfect for on the go practitioners our strong & luxurious portable massage chair requires no tools for assembly and folds up easily.

Features:

- Strongly made steel frame
- Removable chest pad, seat pad, knee support, face rest
- Comfortable design with High density foam
- Luxurious washable & durable PU vinyl
- No tools needed for assembly
- Folds up easily

Specifications:

- Open dimensions: 54x76x120-127cm
- Folded dimensions: 54x76x118cm
- Chair weight: 15kg
- Weight capacity: 180kg
- Warranty: 2 Years on frame

ORDER CODE:
PMCBL - BLACK

[Click For More Information](#)

ULTRASONIC CLEANERS

ULTRASONIC CLEANER 2 LITRE



Our 2L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 2L
- Tank Dimensions: 150 x 140 x 100 (LxWxH)
- External Dimensions: 175 x 165 x 230 (LxWxH)
- Heating power: 100W
- Ultrasonic power: 50W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: No
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD2

[Click For More Information](#)

ULTRASONIC CLEANER 3 LITRE

Our 3L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 3L
- Tank Dimensions: 240 x 140 x 100 (LxWxH)
- External Dimensions: 330 x 165 x 230 (LxWxH)
- Heating power: 100W
- Ultrasonic power: 100W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: No
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD3

[Click For More Information](#)



ULTRASONIC CLEANER 6 LITRE

Our 6L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

Capacity: 6L



Tank Dimensions: 300 x 150 x 150 (LxWxH)
External Dimensions: 400 x 175 x 305 (LxWxH)
Heating power: 300W
Ultrasonic power: 150W
Ultrasonic frequency: 40KHz
Temp Adjustable: Ambient to 80°C
Drainage Valve: Yes
Display: Digital
Timer: 1 - 99min

ORDER CODE:**USD6**[Click For More Information](#)

ULTRASONIC CLEANER 9 LITRE

Our 9L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

Turbo power - normal and soft functions
Auto-degas function
Automatic ultrasonic frequency tracking
Memory function
LED display for temperature and timer
User friendly and clear panel for easy operation
Comes with basket and lid as standard

Specifications:

Capacity: 9L
Tank Dimensions: 300 x 235 x 150 (LxWxH)
External Dimensions: 400 x 265 x 305 (LxWxH)
Heating power: 300W
Ultrasonic power: 200W
Ultrasonic frequency: 40KHz
Temp Adjustable: Ambient to 80°C
Drainage Valve: Yes
Display: Digital
Timer: 1 - 99min

ORDER CODE:**USD9**[Click For More Information](#)



ULTRASONIC CLEANER 13 LITRE

Our 13L ultrasonic cleaners are quiet with variable rinse times, temperatures, stainless steel finish and a clear digital display they are fitting for any modern clinic.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 13L
- Tank Dimensions: 330 x 300 x 150 (LxWxH)
- External Dimensions: 360 x 330 x 310 (LxWxH)
- Heating power: 300W
- Ultrasonic power: 400W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: Yes
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD13

[Click For More Information](#)

ULTRASONIC CLEANER 20 LITRE

Our 20L ultrasonic cleaners are quiet with variable rinse times, temperatures, stainless steel finish and a clear digital display they are fitting for any modern clinic.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:



Capacity: 20L
Tank Dimensions: 500 x 300 x 150 (LxWxH)
External Dimensions: 550 x 330 x 310 (LxWxH)
Heating power: 500W
Ultrasonic power: 400W
Ultrasonic frequency: 40KHz
Temp Adjustable: Ambient to 80°C
Drainage Valve: Yes
Display: Digital
Timer: 1 - 99min

ORDER CODE:**USD20**[Click For More Information](#)

ULTRASONIC CLEANER 27 LITRE



Our 27L ultrasonic cleaners are quiet with variable rinse times, temperatures, stainless steel finish and a clear digital display they are fitting for any modern clinic.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

Capacity: 27L
Tank Dimensions: 500 x 300 x 200 (LxWxH)
External Dimensions: 550 x 330 x 360 (LxWxH)
Heating power: 500W
Ultrasonic power: 500W
Ultrasonic frequency: 40KHz
Temp Adjustable: Ambient to 80°C
Drainage Valve: Yes
Display: Digital
Timer: 1 - 99min

ORDER CODE:**USD27**[Click For More Information](#)

WATER DISTILLERS

WATER DISTILLER BENCHTOP

This compact water distiller provides a suitable water supply solution for many practice autoclaves.

Specifications:

Stainless Steel Tank

750w Power

Produces 1L Per Hour

Internal Dimensions: 180x200mm

External Dimensions: 290x290x390mm

Weight: 3.5KG

CE Mark

ORDER CODE:

WDIS

[Click For More Information](#)

WATER DISTILLER CHARCOAL FILTER

Replacement carbon filter for the WDIS water distiller.

Replacement carbon filter for the WDIS water distiller.

ORDER CODE:

WDISCF

[Click For More Information](#)





WATER DISTILLER CLEANING CRYSTALS

A tub of crystals used in the cleaning of benchtop water distillers.

A tub of crystals used in the cleaning of benchtop water distillers.

Order Code:

WDISCC

[Click For More Information](#)

WATER DISTILLER REPLACEMENT BOTTLE

Replacement water bottle for the WDIS water distiller.

Replacement water bottle for the WDIS water distiller.

Order Code:

WDISB

[Click For More Information](#)



WHEELCHAIRS



STANDARD WHEELCHAIR 18"

A strong economical wheelchair with 110kg capacity. A great addition to any clinic.

Features:

- Single cross brace powder coated frame
- Front rubber tires
- Pneumatic 24" spoke wheel rear tires
- Flip back desk length armrest
- Swing away footrest with plastic footplate

Specifications:

- Total Length 103cm
- Total Width 63cm
- Total Height 91cm
- Width Folded 26.5cm
- Seat Width 45cm
- Seat Depth 40cm
- Backrest Height 44cm
- Seat Height 50cm
- Weight Capacity 110kg
- Front Wheel Diameter 8" x 6"
- Rear Wheel Diameter 24"
- Gross Weight 22.5kg
- Net Weight 20Kg
- Carton Size 82 x 27 x 43cm

ORDER CODE: 18WCS

[Click For More Information](#)

PATIENT MOVER WHEELCHAIR 18"

Our patient mover wheelchair is a great aid to ferry patients around the largest of facilities with ease.

Features:

Single cross brace, powder coating
Front PVC tires
Rear PU tires
Flip back desk length armrest
Swing away footrest with plastic footplate

Specifications:

Total Length 97cm
Total Width 63.5cm
Total Height 92cm
Width Folded 27cm
Seat Width 46cm
Seat Depth 41cm
Backrest Height 43cm
Seat Height 49cm
Weight Capacity 110kg
Front Wheel Diameter 8"
Rear Wheel Diameter 12"
Gross Weight 16.4kg
Net Weight 14.4kg
Carton Size 72 x 25 x 93

ORDER CODE:
18WCM

[Click For More Information](#)



BARIATRIC WHEELCHAIR 22"

Every practice needs a strong wheelchair to carry heavy patients. With a foldable design, 200kg weight capacity and comfortable seat it presents great value.

Features:

Single cross brace powder coated frame
Front PU tires
Spoke rear wheel with PU tires
Flip back desk length armrest
Swing away footrest with plastic footplate
Nylon with foam seat & backrest

Specifications:

Total Length: 103cm
Total Width: 73cm
Total Height: 91cm
Width Folded: 26.5cm
Seat Width: 55cm
Seat Depth: 40cm
Backrest Height: 44cm
Seat Height: 50cm
Weight Capacity: 200kg
Front Wheel Diameter: 8" x 2"
Rear Wheel Diameter: 24"
Gross Weight: 23kg
Net Weight: 20.5kg
Carton Size: 82 x 27 x 93

ORDER CODE:
22WCB

[Click For More Information](#)

THERMOMETRY

INFRARED FOREHEAD THERMOMETER



This infrared forehead thermometer allows for fast, easy and accurate measurement. With a measurement time of one second at a distance of up to 5cm it is suited for use in many clinical settings such as GP, Hospital, Workplace, First Aid and more.

Specifications:

Measurement range: Human body temperature mode: 32.0°C-43.0°C

Display Resolution: 0.1°C

Measurement place: Forehead

Accuracy: 32.0°C - 34.9°C, maximum permissible error $\pm 0.3^{\circ}\text{C}$ / 35.0°C - 42.0°C, maximum permissible error $\pm 0.2^{\circ}\text{C}$ / 42.1°C - 42.9°C, maximum permissible error $\pm 0.3^{\circ}\text{C}$

Display screen: LCD screen with 4-digit and special icons

Operating environment: 10°C-40°C

Max. Relative humidity: <85%

Storage environment: -20°C - +55°C

Battery: 2 x 1.5V AAA Batteries

Order Code:

CYT1

[Click For More Information](#)

BLOOD PRESSURE



BPLAB STANDARD 24 HOUR BLOOD PRESSURE MONITOR

The BPLab® Standard is a compact and easy to use device which provides comprehensive 24 hour blood pressure tests. The device has a confirmed accuracy class A/A for measuring blood pressure according to the BHS-93 protocol in general, paediatric and pregnant women populations. Utilising modern technology of pulse wave assessment this monitor has been designed for the non-invasive measurement of many advanced parameters. The adaptive inflation feature enhances patient comfort with minimum measurement time and Pulse Wave Analysis(PWA) without an additional inflation step.

Medicare Rebate Information:

Item number 11607 will be able to be claimed as of 01/11/2021. See Medicare rebate tab for full information.

MBS Online Link -

http://www9.health.gov.au/mbs/fullDisplay.cfm?type=item&q=11607&q_t=item&criteria=11607

Unique Features:

Measurement Quality Control (Due To The Recorded Oscillograms)
Patient Activity & Body Position Sensor (Optional Extra)
BPStat (Special Software For Tabulation Of Arrayed Data, required for Clinical Studies etc)
Complete Integration Telemedicine Systems Inc Tholomeus (Remote Reporting on 24 Hour ABPMK & PWA)

Specifications:

Measurement Type: Oscillometric
Measured Parameters: Systolic (SYS), Diastolic (DIA), Mean (MBP), Pulse Pressure (PP), Central Aortic Pressure (CAP [SYSao, DIAao, MBPao]), Brachial & Aortic Augmentation Index (AIx), Ejection Duration (ED), Pulse Pressure Amplification (PPA), Maximum Rate Of Blood Pressure Rise (dp/dtmax), Subendocardial Viability Ratio (SEVR) and as a set of options Ankle Brachial Index (ABI) and information on patients activity and body position.

Please note some parameters are only available for view with the purchase of additional software packages.

Warranty: 3 Years

Endorsements:

BHS-93 Protocol A/A grade in general, pediatric and pregnant women population. CAP, AIx: Pass in general population.
Enlisted in BHS, Stride BP & Dabl Educational Trust websites.

Available Cuff Sizes (Latex Free):

Small (16-24CM)
Medium (24-32cm)
Large (32-42cm)
Extra Large (38-50cm)

Standard Package Inclusions:

Basic Package: BPLab® Standard Device, Pouch, Shoulder and Waist Belts, Measuring Tape, Medium Cuff (24-32cm), Large Cuff (32-42cm), Cuff Hose, Cuff Clip, USB Connection Cable, BPLabWin Software Disc

BPLabWin Software: Processing of ABPM data (SYS, DIA, MBP, PP, HR & Double Product)

Order Codes:

BPLABSTD: BPLAB Standard 24 Hour Blood Pressure Monitor

Cuffs:

BPLABPCSMML: Ambulatory Blood Pressure Cuff Small 16-24cm

**BPLABPCMED: Ambulatory Blood Pressure Cuff Medium
24-32cm**

BPLABPCLGE: Ambulatory Blood Pressure Cuff Large 32-42cm

**BPLABPCXLGE: Ambulatory Blood Pressure Cuff Extra Large
38-50cm**

BPLABPCCOV: Cotton Cuff Cover Set Of 5

Accessories & Cables:

**BPLABSTDBAG: Shoulder Carry Bag & Waist Belt for BPLAB
Standard**

BPLABPSTDCAB: PC Connection Cable For BPLAB Standard

Software:

**BPLABVASOTENSOFFICE: Vasotens Office Software (Requires
Software Dongle)**

**BPLABVASOTENS24: Vasotens 24 Software (Requires Software
Dongle)**

[Click For More Information](#)

PULSE OXIMETRY



PULSE OXIMETER FINGERTIP

A small pulse oximeter intended for measuring the pulse oxygen saturation, pulse rate and perfusion index through a users finger. Suitable for monitoring of adult & paediatric patients in medical clinics, hospitals and at home.

Pulse Oximeter Features:

- Special splash proof and drop resistant design
- OLED display of SpO2, PR, PI, pulse bar & waveform
- Spot check and continuous pulse oximeter measurement modes
- Pulse rate analysis for spot check measurement
- Up to 12 groups SpO2 data storage
- Audible & visual alarm with low battery indication
- Automatic power on/off
- Four direction display
- Artefact removal & anti motion
- Accurate & sensitive

Pulse Oximeter Specification:

Spo2

Transducer: Dual-wavelength LED sensor
Measuring range: 35-100%
Measuring accuracy: $\leq 2\%$ range from 70-100%

Pulse Rate

Measuring range: 30-240bpm
Measuring accuracy: ± 2 bpm or $\pm 2\%$ (whichever is the greater)

Perfusion Index

Display range: 0-20%

Default Alarm Limit

SpO2 low limit: 90%
Pulse rate high limit: 120bpm
Pulse rate low limit: 50bpm

Measuring Mode

30 seconds spot check mode
Continuous mode

Dimensions, Weight & Power

Dimensions: 56x34x30mm (L x W x H)
Weight: 52g
Power: 2 x AAA alkaline battery

ORDER CODE:

PMFPO

[Click For More Information](#)

NEW PRODUCTS



NITROSPRAY CRYOTHERAPY FLASK

Premier Nitrospray cryotherapy flask constructed of a double-walled vacuum stainless steel container to hold liquid nitrogen. The outflow track extension port, located on the top of the instrument has a luer-lock connector to easily engage & disengage spray tips and other attachments as required.

Features:

- Double walled vacuum cryotherapy flask
- Easy connect luer lock connector
- Relief valve precisely regulates pressure
- Curved trigger provides unobstructed view of treatment site
- Easy to hold grooved head
- Durable stainless steel construction
- Spray tip protector prevents accidental damage and contact with patient

Specification:

- Holding Time: 6-8 Hours
- Intermittent Use: 3-6 Hours Treatment Time
- Package Inclusions: Flask, Set of five spray tips & Spray tip protector
- Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)
- Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

ORDER CODES:

LNSLITEPLUS - 300ml

LNSPLUS - 475ml

[Click For More Information](#)

CRYOTHERAPY PACKAGE ONE

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle, Extractor & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle, Extractor & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal



Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

Order Code:

LNPACK1

[Click For More Information](#)

CRYOTHERAPY PACKAGE TWO

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal



Dewar Specification:

Capacity: 20L
 Weight (empty): 9.5kg
 Caliber (neck opening of dewar): 50mm
 Outer diameter of dewar: 394mm
 Height of dewar: 672mm
 Static evaporation rate: 0.12L/d
 Static holding time: 168 days
 Ladle outer diameter: 38mm
 Ladle height: 120/276mm
 No. ladles per dewar: 3~6

Order Code:**LNPACK2**[Click For More Information](#)**PULSE OXIMETER FINGERTIP**

A small pulse oximeter intended for measuring the pulse oxygen saturation, pulse rate and perfusion index through a users finger. Suitable for monitoring of adult & paediatric patients in medical clinics, hospitals and at home.

Pulse Oximeter Features:

Special splash proof and drop resistant design
 OLED display of SpO2, PR, PI, pulse bar & waveform
 Spot check and continuous pulse oximeter measurement modes
 Pulse rate analysis for spot check measurement
 Up to 12 groups SpO2 data storage
 Audible & visual alarm with low battery indication
 Automatic power on/off
 Four direction display
 Artefact removal & anti motion
 Accurate & sensitive

Pulse Oximeter Specification:**Spo2**

Transducer: Dual-wavelength LED sensor
 Measuring range: 35-100%
 Measuring accuracy: $\leq 2\%$ range from 70-100%

Pulse Rate

Measuring range: 30-240bpm
 Measuring accuracy: $\pm 2\text{bpm}$ or $\pm 2\%$ (whichever is the greater)

Perfusion Index

Display range: 0-20%

Default Alarm Limit

SpO2 low limit: 90%
 Pulse rate high limit: 120bpm



Pulse rate low limit: 50bpm

Measuring Mode

30 seconds spot check mode
Continuous mode

Dimensions, Weight & Power

Dimensions: 56x34x30mm (L x W x H)
Weight: 52g
Power: 2 x AAA alkaline battery

ORDER CODE:

PMFPO

[Click For More Information](#)

**PML1 LED EXAMINATION LIGHT**

The PML1 LED examination light is a great all-round diagnostic light. It provides a soft white light for better colour rendition and easier diagnosis and treatment. With a small light head and convenient handle with integrated on/off switch it is easy to move without blocking the vision of the operator. The flexible gooseneck allows positioning in a diverse range of angles for optimal illumination of the work area. A stepless dimmer control allows the operator to set the precise illumination required for each use without being constrained to factory preset levels and the heat free LED provides bright white light for up to 50,000 hours of operation.

It is designed for use in Consulting Rooms, Recovery Areas, Emergency Rooms, Gynaecology Examinations and Minor Surgical Procedures.

Features:

- Powerful long-life LED with up to 50,000 hours of operation
- Low energy consumption
- IR and heat free cold white light
- Comfortable guiding handle with integrated On/Off switch
- Cooling through a heatsink so no fans or moving parts
- Superb colour rendering
- Mobile base and wall mount provided with every light

Specifications:

- Light Source: 3W LED
- Light Intensity: 89,000 Lux @ 30cm / 40,000 Lux @ 50cm
- Light Field Diameter: 80-150mm
- Colour Temperature: 4,500-5,000K
- Lamp Head Diameter: 75mm
- Lamp Head Length: 110mm
- Base Diameter: 530mm
- Base Height: 790mm
- Gooseneck Length: 760mm
- Total Height: 1,830mm

Order Codes:

PML1WH - WHITE**PML1BL - BLACK**[Click For More Information](#)

PML2 LED EXAMINATION LIGHT

The PML2 LED examination light has 6 powerful LED's offering powerful illumination and great value. It provides white light for better colour rendition and easier diagnosis and treatment. With a small light head with integrated handle it is easy to move without blocking the vision of the operator. The flexible goose neck allows positioning in a diverse range of angles for optimal illumination of the work area. A stepless dimmer control allows the operator to set the precise illumination required for each use without being constrained to factory preset levels and the heat free LED provides bright white light for up to 50,000 hours of operation.

It is designed for use in Consulting Rooms, Recovery Areas, Emergency Rooms, Gynaecology Examinations & Minor Surgical Procedures.

Features:

- Powerful long-life LED with up to 50,000 hours of operation
- Low energy consumption
- IR and heat free cold white light
- Small light head with integrated handle
- Cooling through a heatsink so no fans or moving parts
- Superb colour rendering
- Mobile base and wall mount provided with every light

Specifications:

- Light Source: 6W LED
- Light Intensity: 78,000 Lux @ 30cm / 40,000 Lux @ 50cm / 20,000 Lux @ 100cm
- Light Field Diameter: 70-200mm
- Colour Temperature: 4,500-5,500K
- Lamp Head Diameter: 80mm
- Base Diameter: 530mm
- Base Height: 790mm
- Gooseneck: Length 760mm
- Total Height: 1,815mm

Order Codes:

PML2WH - WHITE**PML2BL - BLACK**[Click For More Information](#)



BPLAB STANDARD 24 HOUR BLOOD PRESSURE MONITOR

The BPLab® Standard is a compact and easy to use device which provides comprehensive 24 hour blood pressure tests. The device has a confirmed accuracy class A/A for measuring blood pressure according to the BHS-93 protocol in general, paediatric and pregnant women populations. Utilising modern technology of pulse wave assessment this monitor has been designed for the non-invasive measurement of many advanced parameters. The adaptive inflation feature enhances patient comfort with minimum measurement time and Pulse Wave Analysis(PWA) without an additional inflation step.

Medicare Rebate Information:

Item number 11607 will be able to be claimed as of 01/11/2021. See Medicare rebate tab for full information.

MBS Online Link -

http://www9.health.gov.au/mbs/fullDisplay.cfm?type=item&q=11607&q_t=item&criteria=11607

Unique Features:

Measurement Quality Control (Due To The Recorded Oscillograms)
 Patient Activity & Body Position Sensor (Optional Extra)
 BPStat (Special Software For Tabulation Of Arrayed Data, required for Clinical Studies etc)
 Complete Integration Telemedicine Systems Inc Tholomeus (Remote Reporting on 24 Hour ABPMK & PWA)

Specifications:

Measurement Type: Oscillometric
 Measured Parameters: Systolic (SYS), Diastolic (DIA), Mean (MBP), Pulse Pressure (PP), Central Aortic Pressure (CAP [SYSao, DIAao, MBPao]), Brachial & Aortic Augmentation Index (AIx), Ejection Duration (ED), Pulse Pressure Amplification (PPA), Maximum Rate Of Blood Pressure Rise (dp/dtmax), Subendocardial Viability Ratio (SEVR) and as a set of options Ankle Brachial Index (ABI) and information on patients activity and body position.

Please note some parameters are only available for view with the purchase of additional software packages.

Warranty: 3 Years

Endorsements:

BHS-93 Protocol A/A grade in general, pediatric and pregnant women population. CAP, AIx: Pass in general population.
 Enlisted in BHS, Stride BP & Dabl Educational Trust websites.

Available Cuff Sizes (Latex Free):

Small (16-24CM)
 Medium (24-32cm)
 Large (32-42cm)
 Extra Large (38-50cm)

Standard Package Inclusions:

Basic Package: BPLab® Standard Device, Pouch, Shoulder and Waist

Belts, Measuring Tape, Medium Cuff (24-32cm), Large Cuff (32-42cm), Cuff Hose, Cuff Clip, USB Connection Cable, BPLabWin Software Disc

BPLabWin Software: Processing of ABPM data (SYS, DIA, MBP, PP, HR & Double Product)

Order Codes:

BPLABSTD: BPLAB Standard 24 Hour Blood Pressure Monitor

Cuffs:

BPLABPCSMML: Ambulatory Blood Pressure Cuff Small 16-24cm

BPLABPCMED: Ambulatory Blood Pressure Cuff Medium 24-32cm

BPLABPCLGE: Ambulatory Blood Pressure Cuff Large 32-42cm

BPLABPCXLGE: Ambulatory Blood Pressure Cuff Extra Large 38-50cm

BPLABPCCOV: Cotton Cuff Cover Set Of 5

Accessories & Cables:

BPLABSTDBAG: Shoulder Carry Bag & Waist Belt for BPLAB Standard

BPLABPSTDCAB: PC Connection Cable For BPLAB Standard

Software:

BPLABVASOTENSOFFICE: Vasotens Office Software (Requires Software Dongle)

BPLABVASOTENS24: Vasotens 24 Software (Requires Software Dongle)

[Click For More Information](#)

INDEX

1

10L Liquid Nitrogen Dewar 81

2

20L Liquid Nitrogen Dewar 82

3

30L Liquid Nitrogen Dewar 83

A

Anaesthesia Trolley 60

B

Bariatric Wheelchair 22" 103

Beauty Day Spa Massage Couch 18

Blood Collection Chair 89

BPLab Standard 24 Hour Blood Pressure Monitor 106, 115

C

Chiropractic Premium Fixed Height Table 11

Chiropractic Premium Table 10

Chiropractic Table 9

Contoured Massage Couch 15

Control Box For Electric Couches 29

Cryotherapy Package One 77, 84, 111

Cryotherapy Package Two 77, 84, 111

D

Difficult Airway Trolley 63

Double Linen Skip 49

Double Linen Skip With Lid 51

Double Step Stool 92

Dressing Change Trolley 60

E

Emergency Trolley 56

Extractor for 10L Dewar 86

Extractor For 20L & 30L Dewars 86

F

Fitted Disposable Bedsheets BX100 24

Five Section Treatment Couch No Postural Drainage 12

Five Section Treatment Couch With Postural Drainage 13

FlexLED Examination Light 68

Foot Control Red Dot For Electric Couch Hi-Lo 28

Foot Control Yellow Dot For Electric Couch Tilt 28

Four Section Physio Traction Treatment Couch 17

Free Standing Treatment Couch 21

Free Standing Treatment Couch With Drawers 22

G

Gas Strut Long For Electric Couches 27

Gas Strut Short For Electric Couch 26

Gynaecology Premium Treatment Couch 8
Gynaecology Treatment Couch 7

H

Hand Control For Electric Couch Hi-Lo 27

I

Infrared Forehead Thermometer 105
Infusion Trolley 61
Instrument Trolley No Shelf 62
Instrument Trolley One Shelf 62

L

Ladle For 10L Dewars 87
Ladle For 20L & 30L Dewars 87
Linen Skip Bag Soluble Seam BX250 52
Linen Skip Bag Standard 52
Liquid Nitrogen Gloves 88

M

Mayo Instrument Trolley 47
Medical Computer Trolley With Printer Shelf 65
Medical Records Trolley No Drawer 63
Medical Records Trolley One Drawer 64
Medicine Trolley 57

N

Neurological Bobath Treatment Couch 14
Nitrospray Accessory Adaptor For Screw Tips 79
Nitrospray Cryotherapy Flask 76, 110
Nitrospray Cryotherapy Spray Tips 78
Nitrospray Needle Extender Rigid 80
Nitrospray Spray Tip Extension Bent 80
Nitrospray Spray Tip Protector 79
Non-Fitted Disposable Bedsheets BX100 24

O

Oxygen Trolley Size D 48

P

Paper Towel Holder For Electric & Fixed Height Couches 26
Patient Mover Wheelchair 18" 103
Phantom LED Minor Procedure Light 72
PML1 LED Examination Light 66, 113
PML2 LED Examination Light 67, 114
Podiatry / Multi Purpose Chair 6
Portable Massage Chair 94
Portable Massage Couch 23
Powder Coated Trolley Four Drawer 55
Powder Coated Trolley No Drawer 53
Powder Coated Trolley One Drawer 53
Powder Coated Trolley Three Drawer 54
Powder Coated Trolley Two Drawer 54
Pulse Oximeter Fingertip 108, 112

R

Round Stool Premium With Foot Ring 90
Round Stool Standard 90

S

Saddle Stool Premium With Foot Ring 91
 Saddle Stool Standard 91
 Single Linen Skip 49
 Single Linen Skip With Lid 50
 Single Section Treatment / Change Couch With Side Rails 20
 Single Step Stool 92
 Slimline X-Ray Viewers 74
 Splitter For Electric Couches 29
 Stainless Steel Anaesthetic Trolley Four Drawer 43
 Stainless Steel Anaesthetic Trolley Six Drawer 44
 Stainless Steel Double Trolley Four Drawer (Two Over Two) 39
 Stainless Steel Double Trolley No Drawer 37
 Stainless Steel Double Trolley One Drawer (Full Width) 37
 Stainless Steel Double Trolley Two Drawer (Full Width) 38
 Stainless Steel Double Trolley Two Drawer (Side By Side) 39
 Stainless Steel IV Stand Four Hook 46
 Stainless Steel IV Stand Two Hook 46
 Stainless Steel Resuscitation Trolley Five Drawer 44
 Stainless Steel Resuscitation Trolley Seven Drawer 45
 Stainless Steel Trolley Four Drawer 33
 Stainless Steel Trolley Four Drawer With Top Locking Drawer 35
 Stainless Steel Trolley No Drawer 30
 Stainless Steel Trolley No Rails 50x50x90 41
 Stainless Steel Trolley No Rails 60x50x90 41
 Stainless Steel Trolley No Rails 80x50x90 42
 Stainless Steel Trolley One Drawer 31
 Stainless Steel Trolley One Drawer With Lock 33
 Stainless Steel Trolley Three Drawer 32
 Stainless Steel Trolley Three Drawer With Top Locking Drawer 35
 Stainless Steel Trolley Two Drawer 31
 Stainless Steel Trolley Two Drawer With Top Locking Drawer 34
 Standard Wheelchair 18" 102
 Standard X-Ray Viewers 75

T

Three Section Bariatric Treatment Couch 3
 Three Section Medical Treatment Couch 2
 Three Section Medical Treatment Couch With Side Rails 19
 Three Section Physio Treatment Couch 5
 Tilt Table Two Section 16
 Treatment Trolley One Drawer 58
 Treatment Trolley Three Drawer 59
 Treatment Trolley Two Drawer 58
 Triple Linen Skip 50
 Triple Linen Skip With Lid 51
 Trolley For 20L & 30L Dewars 88
 Two Section Treatment Couch 4

U

ULED Examination Light 69
 ULED Plus Examination Light 70
 Ultrasonic Cleaner 13 Litre 98
 Ultrasonic Cleaner 2 Litre 95
 Ultrasonic Cleaner 20 Litre 98
 Ultrasonic Cleaner 27 Litre 99
 Ultrasonic Cleaner 3 Litre 96

Ultrasonic Cleaner 6 Litre 96

Ultrasonic Cleaner 9 Litre 97

W

Water Distiller Benchtop 100

Water Distiller Charcoal Filter 100

Water Distiller Cleaning Crystals 101

Water Distiller Replacement Bottle 101



Maquarie Medical Systems

E sales@machealth.com.au

P (02) 9692 7921

www.macquariemed.com.au